



Dermal

Podiatry

Instruments



800.221.2166 | www.sklarcorp.com



© Copyright 2016.
Sklar Corporation, West Chester, PA. All rights reserved. No portion of this catalog may be reproduced without written permission of Sklar.

TABLE OF CONTENTS

Nail Nippers & Splitters	3-7
Bone Forceps & Rongeurs	8-12
Curettes.....	13-16
Extractors	17
Punches.....	18
Bone Files & Rasps	19-21
Elevators	22-24
Osteotomes	25
Hand Instruments	26-28
Forceps	29-52
Scissors	53-77
Needle Holders	78-84
Retractors	85-90
Pin & Wire	91-94
Plaster Instruments	95-97
Diagnostics	98-99
Blade Breakers & Ring Cutters.	100
Disposables	101-106
Instrument Sets	107-114
Care & Cleaning.....	115-129



DISCLAIMERS:

Sklar offers a full line of quality surgical products and services to meet your needs.

Product information within this catalog is intended to provide general information about Sklar, its products and services and makes no representation either expressed or implied to the contrary. It is not meant to provide specialist advice or instructions concerning the products available from Sklar. For additional information and special inquiries concerning our products and services, please contact Sklar directly.

In an effort to continually improve our product line, Sklar reserves the right to change patterns without notice. Contact your sales representative to verify availability.

In no event is Sklar liable for any incidental, special or exemplary damages or loss of profits in connection with the use of a Sklar Instrument or Product.

Convenience Kits (FDA recognized term) illustrated in this catalog are representative of Sklar's assembly capability. The final design and approval of the kit is the responsibility of the customer and will be managed in accordance with the Sklar Design and Development process.

Sklar Corporation
889 S. Matlack Street
West Chester, PA 19382

Phone: 800.221.2166
Fax: 610.696.9007
Email: Surgi@sklarcorp.com

ORDERS

- Sklar's full line of Instrument Products are available through any authorized Sklar Distributor. Contact our Customer Service Department at 800.221.2166 to locate a Sales Representative in your area.
- Customers with established accounts and in good standing are subject to Net 30 days from date of invoice. Poor credit conditions can result in denial, delay or cancellation of orders.

SHIPMENTS

- Pricing is FOB West Chester PA unless otherwise noted and may be subject to change. Shipping schedules are approximate and subject to change without notice. Same day shipping is available. Standard orders will be shipped via the carrier of Sklar's choice. Special shipping requests will be accepted but additional handling charges are applicable.
- All Sklar instruments and product are carefully inspected prior to shipment. Customers shall inspect all items upon receipt and file a claim within (10) days of receipt of shipment.


WARRANTY

- Instruments and Products are recommended for their intended use only, and require continuous and proper care. Misuse can damage instruments and products and void the warranty.
- All Sklar instruments are guaranteed to be free from defect in material and workmanship. Defective instruments may be returned for repair or replacement at no charge. Repair or replacement is at the discretion of Sklar. Warranty is void on any instrument serviced by any company other than Sklar Instrument Repair Services.
- Sklar offers a variety of supplementary surgical products. These products are held to high standards for material and craftsmanship and some are subject to a limited warranty (vary by product).

RETURNS

- Unused instruments returned in the original packaging will be accepted within 30 days of purchase. Standard sterile kits may be returned within 30 days provided sterility of the inner packs have not been compromised, package is unopened and resalable. Most cleaning products are non-returnable.
- Product that was special (engraved/etched) or custom ordered, discontinued, marked or altered in any way, including product or package labeling, will not be accepted for return or credit.
- A Returned Goods Authorization (RGA) number must be obtained prior to returning products. Contact Customer Service at 800.221.2166 for authorization. Credit, repair and/or replacement are at Sklar's discretion. See www.sklarcorp.com for complete return policy and procedures.






LITTAUER CUTTING FORCEPS

STRAIGHT BACK-STAINLESS


97-1251	4"
97-1255	4½"
97-1259	5"
97-1262	6"
97-1261	6¼"



SKLARLITE™ LITTAUER CUTTING FORCEPS

STRAIGHT BACK-STAINLESS


23-9251	4"
23-9255	4½"
23-9259	5"



LITTAUER CUTTING FORCEPS

CURVED BACK-STAINLESS

97-1280	4½"
97-1285	5"
97-1286	6"



SKLARLITE™ LITTAUER CUTTING FORCEPS

CURVED BACK-STAINLESS


23-9273	4½"
23-9275	5"
23-9276	6"



SKLAR-TAPER™ LITTAUER CUTTING FORCEPS

STAINLESS


97-1174 5"



SKLAR LITTAUER CUTTING FORCEPS

DOUBLE-ACTION


97-1185 6" Straight



NAIL NIPPERS

CONCAVE BLADES


97-1308 6" Stainless Narrow
97-1309 5½" Stainless Heavy



NAIL NIPPERS

CONCAVE BLADES- STAINLESS


97-1145 4½" Folding Spring
97-1148 4½" Double Spring




NAIL NIPPER

CONCAVE BLADES


97-1055 5½" Stainless




NAIL NIPPER
97-1056 5½" Straight Back




NAIL NIPPERS
CONCAVE BLADES-BARREL SPRING
97-1045 4½"
97-1046 5½"



MYCOTIC NAIL NIPPER
DOUBLE ACTION
97-1057 5¾" Stainless



NAIL NIPPER
ANGLED CONCAVE BLADES
97-1156 5½" Stainless



TISSUE NIPPERS

CONVEX BLADES – STAINLESS


97-1241	4"	6mm
97-0335	4"	9mm
97-1242	4½"	11mm
97-1243	5"	14mm



SKLARLITE™ TISSUE NIPPER

CONVEX BLADES–STAINLESS

23-9292	4½"	Double Spring	11mm
---------	-----	---------------	------



TISSUE FORCEPS

SCISSORS HANDLE–STRAIGHT JAWS


97-1238	Stainless	4"
---------	-----------	----



CUTICLE AND NAIL NIPPERS

RING HANDLE


97-1235	3½"	Finger Nail	Straight
97-1245	4¼"	Toe Nail	Curved



NAIL SPLITTERS

TAPERED JAW, STRAIGHT BACK-STAINLESS


97-1660	4"
97-1661	4½"
97-1662	5"
97-1663	6"



PLATYPUS NAIL FORCEPS

PLIER HANDLE


97-1318	5½"	6mm wide
97-1315	5½"	4mm wide



INGROWN NAIL SPLITTING SCISSOR

ONE POINTED BLADE


97-0555	6"	One Blade Serrated
---------	----	--------------------



INGROWN TOE NAIL FORCEPS

ENGLISH (ANVIL) PATTERN


97-1150	5"	Stainless
---------	----	-----------



SKLARLITE™ INGROWN TOE NAIL FORCEPS


ENGLISH (ANVIL) PATTERN

23-9250	5"	Stainless
---------	----	-----------



PNA FORCEPS

97-1340 4½"



LISTON BONE FORCEPS

SHARP POINTS


40-3040 5½" Straight



ROWLAND BONE FORCEPS

DOUBLE ACTION


41-1175 7" Narrow Jaws



MCINDOE BONE CUTTING FORCEPS

DOUBLE ACTION


41-1170 7½" Slight Angle



LISTON BONE FORCEPS

STRAIGHT


40-3055	5½"
40-3077	7¾"
40-3087	8¾"



LISTON BONE FORCEPS

ANGULAR

40-3155	5½"
40-3167	6¾"
40-3177	7¾"
40-3187	8¾"



KLEINERT-KUTZ BONE CUTTING FORCEPS

DOUBLE ACTION

40-4560	6"	Straight
40-4562	6"	Delicate



RUSKIN-LISTON BONE CUTTING FORCEPS


DOUBLE ACTION

40-4572	7¼"	Straight
40-4575	7¼"	Up-Angled




RONGEUR NO. 5S

40-3196 5 3/4" Side Cut



RONGEUR NO. 5

40-3198 6 1/2" Side Cut



MINI BLUMENTHAL RONGEURS


DELICATE

67-5987 4 3/4" 30° Angle
67-5988 4 3/4" 45° Angle




BLUMENTHAL RONGEURS

67-5990 6" 30° Angle
67-5991 7" 30° Angle
67-5992 6" 45° Angle
67-5994 6" 90° Angle



LEMPERT RONGEURS


47-2467 6 1/4" Straight
47-2468 6 1/4" Curved
47-2465 7 1/2" Straight
47-2466 7 1/2" Curved



MICRO FRIEDMAN RONGEURS

1.3MM JAWS, VERY DELICATE


40-3186	5½"	Straight
40-3189	5½"	Curved



FRIEDMAN BONE RONGEURS

DELICATE


40-3190	5½"	Curved
40-3192	6½"	Curved



CLEVELAND RONGEURS


CURVED

40-3200	5½"	#4
40-3202	6½"	#4A




MEAD RONGEURS

40-3206	6½"	Slight Curve
40-3208	6½"	Curved




BANE RONGEUR

67-5175	7"	Curved	3x16mm Bite
---------	----	--------	-------------



HARTMANN RONGEUR


67-5458 5¾" Full Curve 4x16mm Bite



KLEINERT-KUTZ RONGEURS


DOUBLE ACTION

40-4190 6" Slight Curve 2x13mm Bite
40-4192 6" Curved 3x16mm Bite
40-4194 6" Full Curve 3x8mm Bite




ZAUFEL-JANSEN RONGEUR

40-4265 9" 5x15mm Curved Bite



RUSKIN RONGEURS

40-4360 6" Straight 2mm Bite
40-4365 6" Curved 2mm Bite
40-4272 7¼" Straight 5mm Bite
40-4372 7¼" Curved 5mm Bite
40-4472 7¼" Angled to Side 5mm Bite
40-4512 9½" Curved 5mm Bite



CURETTE-EXCAVATORS

DOUBLE END WITH HOLE 5 1/2"


97-0532 #58-1-2 1.5-2mm
97-0536 #58-1-3 1.5-2.5mm



CURETTE-EXCAVATORS

SINGLE END WITH HOLE 5"


97-0512 #58-0 1mm
97-0522 #58-1 1.5mm
97-0524 #58-2 2.0mm
97-0528 #58-3 2.5mm



CURETTE-EXCAVATORS

SINGLE END WITHOUT HOLE 5"


97-0503 #57-0 1mm
97-0505 #57-1 1.5mm
97-0507 #57-2 2mm
97-0509 #57-3 2.5mm



VERUCCA CURETTE

WITHOUT HOLE


97-0564 5" 4mm Cup



VOLKMAN CURETTES

DOUBLE END

40-7250	5"	Oval/Round
40-7265	6 1/2"	Oval/Round
40-7280	8"	Oval/Round



BRUNS CURETTES


STRAIGHT, OVAL

6 3/4"

40-7118	#5-0	40-7112	#2
40-7120	#4-0	40-7113	#3
40-7000	#3-0	40-7114	#4
40-7100	#2-0	40-7115	#5
40-7110	#1-0	40-7116	#6
40-7111	#1		

9"

40-7199	#6-0	40-7205	#1
40-7200	#5-0	40-7206	#2
40-7201	#4-0	40-7207	#3
40-7202	#3-0	40-7208	#4
40-7203	#2-0	40-7209	#5
40-7204	#0	40-7210	#6



BRUNS CURETTES


ANGLED, OVAL

6 1/2"

41-1710	#5-0	41-1716	#2
41-1711	#4-0	41-1717	#3
41-1712	#3-0	41-1718	#4
41-1713	#00	41-1719	#5
41-1714	#0	41-1720	#6
41-1715	#1		

9"


40-7219	#6-0	40-7225	#1
40-7220	#5-0	40-7226	#2
40-7221	#4-0	40-7227	#3
40-7222	#3-0	40-7228	#4
40-7223	#2-0	40-7229	#5
40-7224	#0	40-7230	#6



SKLAR DERMAL CURETTES

OVAL 6 3/4"

06-4112	#0
06-4113	#1
06-4114	#2
06-4115	#3




FOX DERMAL CURETTES

OVAL 5 1/2"

06-4133	3mm
06-4134	4mm
06-4135	5mm
06-4136	6mm

ROUND 5 1/2"


06-4138	3mm
06-4139	4mm
06-4140	5mm
06-4141	6mm



PIFFARD CURETTES


OVAL 6"

06-4130	#1
06-4131	#2
06-4132	#3




REIN CURETTES

OVAL		ROUND	
06-4053	#1	06-4059	2mm
06-4058	#6	06-4060	3mm
		06-4061	4mm




CANNON CURETTES

06-4001	#00
06-4002	#0
06-4003	#1
06-4004	#2
06-4005	#3
06-4006	#4
06-4007	#5




BUCK EAR CURETTES

BLUNT			
STRAIGHT		ANGLED	
67-2488	#000	67-2508	#000
67-2489	#00	67-2509	#00
67-2490	#0	67-2510	#0
67-2491	#1	67-2511	#1
67-2492	#2	67-2512	#2
67-2493	#3	67-2513	#3




BUCK EAR CURETTES

SHARP			
STRAIGHT		ANGLED	
67-2519	#000	67-2498	#000
67-2514	#00	67-2499	#00
67-2515	#0	67-2500	#0
67-2516	#1	67-2501	#1
67-2517	#2	67-2502	#2
67-2518	#3	67-2503	#3



MEYHOEFER CHALAZION CURETTES

65-3169	#00	1mm
65-3170	#0	1.5mm
65-3171	#1	1.8mm
65-3172	#2	2.2mm
65-3173	#3	2.5mm
65-3174	#4	3mm



HEATH CHALAZION CURETTES

ONE SHARP EDGE AND ONE DULL EDGE		
65-3180	#0	.5mm
65-3181	#1	1mm
65-3182	#2	2mm
65-3183	#3	3mm



SKEELE CHALAZION CURETTE


65-3162 2mm Serrated



HEBRA CHALAZION CURETTE

OVAL

65-3185 1mm
65-3186 2mm
65-3187 2.5mm
65-3188 3mm
65-3189 3.5mm



MARTINI BONE CURETTE

DOUBLE END, ROUND CUPS, 5½"


40-7282 2.5x3.5mm
40-7283 3.5x4.5mm
40-7284 4x5mm



WILLIGER OVAL CURETTE


DOUBLE END, OVAL CUPS, 5½"

40-7285 #00/0, 2x2.5mm
40-7287 #0/1, 3x4mm
40-7286 #1/2, 3x4mm




SCHAMBERG COMEDONE EXTRACTORS

- 06-4120 Regular Size (A)
- 06-4121 One End Extra Fine (B)
- 06-4122 Square Ends (C)



UNNA COMEDONE EXTRACTORS

- DOUBLE END**
- 06-4125 5¾" Oval
- 06-4191 5¾" Round



COMEDONE EXTRACTOR SPEAR POINT

- 06-4126 4"




COMEDONE EXTRACTOR SPEAR POINT WITH COVER

- 06-4128 Double End




COMEDONE EXTRACTOR DOUBLE END

- 06-4123 4½" Flat and Wire End Expressors



KEYES CUTANEOUS PUNCHES

- 4"
- | | |
|-------------|-------------|
| 06-4142 2mm | 06-4146 6mm |
| 06-4143 3mm | 06-4147 7mm |
| 06-4144 4mm | 06-4148 8mm |
| 06-4145 5mm | |



KEYES DERMAL PUNCH

- | | |
|---------|--|
| 06-4100 | Set of 6, 2mm thru 7mm with Universal Handle and Sterilizable Base |
| 06-4110 | KEYES Sharpening Stone |

SKLAR TRU-PUNCH™
DISPOSABLE BIOPSY PUNCH

The SKLAR Tru-punch has a seamless steel tip carefully honed for a razor sharp incision and perfect sample. The ribbed and hollow handle permits tissue to be easily dislodged.

STERILE DISPOSABLE BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION	QUANTITIES	PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-0192	1mm	25/bx	96-1092	1mm	50/bx
96-1097	1.5mm	25/bx	96-1098	1.5mm	50/bx
96-1140	2mm	25/bx	96-1100	2mm	50/bx
96-0102	2.5mm	25/bx	96-1102	2.5mm	50/bx
96-1142	3mm	25/bx	96-1105	3mm	50/bx
96-1144	3.5mm	25/bx	96-1110	3.5mm	50/bx
96-1146	4mm	25/bx	96-1115	4mm	50/bx
96-1148	5mm	25/bx	96-1120	5mm	50/bx
96-1150	6mm	25/bx	96-1125	6mm	50/bx
96-1152	8mm	25/bx	96-1130	8mm	50/bx
96-1160	Assorted - 5 of each: 2, 3, 3.5, 4 & 5mm	25/bx	96-0160	Assorted - 5 of each: 2, 3, 3.5, 4 & 5mm	50/bx





NAIL AND BONE RASPS


DOUBLE END

97-1610 Inside/Outside Cut #92

SINGLE END

97-1600 Inside Cut #92A


97-1605 Outside Cut #92B



BONE FILE #12

DOUBLE END


97-1620 7" Plain Serrations



BONE FILE #12A

DOUBLE END

97-1625 7" Plain Serrations & Cross Serrations




NAIL AND BONE RASPS

97-1640 7" Plain Serrations (A)

97-1642 6½" Polokoff Medium Serrations (B)


Double End,
3 mm & 4 mm
Single End



NAIL AND BONE RASP

DOUBLE END

97-1645 7½" Plain Serrations



BONE RASP #10

DOUBLE END—CURVED & FLAT ENDS


97-1615 7½"



BONE FILE

BAYONET HANDLE


40-9685 9"



TC MALTZ RASPS

SINGLE END


97-1505 6¾" Fine
97-1506 6¾" Medium



FOMON RASPS


DOUBLE ENDED, FLAT & CONVEX

40-9686 8¼" 2 Sided, Coarse
40-9688 8½" 4 Sided, Coarse & Fine



LEWIS RASPS


70-2292 7" Fine
70-2293 7" Coarse



JOSEPH RASP


STRAIGHT HANDLE

40-9825 6¼" Cross Serrations




MALTZ RASP

40-9815 7" Backward Cutting



COTTLE RASP

40-9820 8" Backward Cutting



AUFRICHT GLABELLAR RASPS


40-9840 8" Forward Cutting
40-9841 8" Backward Cutting



GALLAHER RASPS

COARSE

70-2295 6½" Straight
70-2296 6½" Curved Left
70-2297 6½" Curved Right



MINI HOHMANN BONE ELEVATORS

40-7450 6½"x6mm
40-7451 6½"x8mm




HOHMANN BONE ELEVATORS

40-7447	8½"	8mm
40-7448	9"	10mm
40-7452	9½"	18mm
40-7453	9½"	43mm
40-7454	10"	24mm
40-7455	10"	20mm
40-7456	10"	22mm
40-7458	10½"	65mm




ELEVATOR

SEMI-SHARP
97-0048 5½" Curved Edge



FROST ELEVATORS

97-0640 Narrow
97-0641 Wide



LOCKE ELEVATORS

97-0623 Narrow 3x15mm
97-0634 Wide 5x20mm



FREER ELEVATOR

DOUBLE END
47-2380 7¼" Sharp/Blunt



SAYRE ELEVATORS

DOUBLE END

40-9265	6½"	Sharp
40-9365	6½"	Blunt



OHL ELEVATOR

DOUBLE END

97-0050	9"
---------	----



McGLAMRY ELEVATORS

CURVED, 6½"


40-7478	9mm
40-7480	11mm
40-7482	13mm
40-7484	15mm
40-7486	17mm



COTTLE PERIOSTEAL ELEVATOR

SLIGHT CURVE

41-1250	7½"	9mm
---------	-----	-----



JOSEPH PERIOSTEAL ELEVATORS


41-1245	6¼"	4.5mm
41-1240	7"	6.5mm



MCKENTY SEPTUM ELEVATORS

CURVED


47-2385	6"	4mm
47-2388	6"	5mm



CUSHING ELEVATORS

SLIGHT CURVE

40-8323	6"	6mm	Hollow Handle
40-8324	7½"	6mm	Solid Handle



LANGENBECK ELEVATORS

SINGLE END

40-8675	7½"	Narrow
40-8775	7½"	Broad



CREGO ELEVATOR

FULL CURVE

40-7410	8"	14mm
---------	----	------



CREGO ELEVATOR

FULL CURVE

40-7415	8"	7mm	Round End
---------	----	-----	-----------



CREGO ELEVATOR

SLIGHT CURVE

40-7420	8"	7mm
---------	----	-----



SELDIN ELEVATOR

97-0527	8"	Double End
---------	----	------------



LAMBOTTE OSTEOTOMES

STRAIGHT, DELICATE, 5"

40-6532 2mm	40-6540 10mm
40-6533 3mm	40-6542 12mm
40-6534 4mm	40-6544 14mm
40-6535 5mm	40-6545 15mm
40-6536 6mm	40-6548 18mm
40-6537 7mm	40-6549 19mm
40-6538 8mm	40-6550 20mm
40-6552 Set of 9, Includes: 2, 4, 6, 8, 10, 12, 14, 18 and 20mm	



JR. LAMBOTTE OSTEOTOMES

STRAIGHT 7"


40-6624 4mm $\frac{5}{32}$ "	40-6632 12mm $\frac{15}{32}$ "
40-6625 5mm $\frac{13}{64}$ "	40-6636 16mm $\frac{5}{8}$ "
40-6626 6mm $\frac{1}{4}$ "	40-6639 19mm $\frac{3}{4}$ "
40-6628 8mm $\frac{5}{16}$ "	40-6642 22mm $\frac{7}{8}$ "
40-6630 10mm $\frac{3}{8}$ "	40-6645 25mm 1"



HOKE OSTEOTOMES

STRAIGHT 5½"


40-6953 3mm	40-6970 10mm
40-6954 4mm	40-6973 13mm
40-6956 6mm	40-6976 16mm
40-6958 8mm	40-6979 19mm



PROBE AND PACKER

DOUBLE END

97-0539 6" Blunt



SPATULA AND PACKER

DOUBLE END

97-0542 5¾"



SPATULA

SQUARE END


97-0580 11-S 6mm



SPATULA

ROUNDED END

97-0581 4mm



EXCAVATOR & PACKER

DOUBLE END


97-0562 6"



PROBE & SPATULA

DOUBLE END


97-0548 6"




EXCAVATOR & PACKER

DOUBLE END


97-0546 5¾"




**INGROWN NAIL SHAVER
FENESTRATED**
97-0554




**BAKER INGROWN
NAIL CUTTER**
97-0552 5"




**TISSUE TWIST
HOOK**
47-2670 5½"




**SPATULA
ROUNDED END**
97-0582 6mm #24A




**SPATULA
ROUNDED END**
97-0583 7mm #24



**NAIL SPLITTERS
STRAIGHT**
97-0601 Narrow 2mm
97-0604 Wide 3.5mm



NAIL SPLITTER
97-0606 Angled, Narrow 2mm




**NUCLEUS KNIFE
SINGLE END**
97-0568 Narrow 2mm




GOLF STICK KNIVES

97-0026	Medium	#2
97-0027	Large	#3



GOLF STICK KNIVES

97-0028	Style A
97-0030	Style B



SOFT CORN KNIVES

ROUND BLADE

97-0035	Straight	4¼"
97-0036	Curved Blade	5"



CHISEL KNIFE

STRAIGHT EDGE

97-0010	Small	8mm
---------	-------	-----



CORN CHISEL

CURVED EDGE

97-0545	5¼"	Spoon Shape
---------	-----	-------------



TC ALLIS TISSUE FORCEPS

ATRAUMATIC

97-1555	6"
97-1560	7½"



TC BABCOCK FORCEPS

ATRAUMATIC

97-1565	6¼"
97-1570	8"



TC ADSON DRESSING FORCEPS

SERRATED


51-3447	4¾"
51-3448	6"



TC ADSON TISSUE FORCEPS

1X2 TEETH, WITH TYING PLATFORM


51-3049	4¾"	Serrated
51-3050	4¾"	Smooth, Gorney-Adson
51-3051	6"	Serrated



TC ADSON-BROWN FORCEPS

SIDE GRASPING

51-1446	4¾"	4x4 Teeth
51-1447	4¾"	7x7 Teeth



TC LITEGRIP™ ADSON-BROWN FORCEPS

SIDE GRASPING

51-1448	4¾"	7x7 Teeth
---------	-----	-----------



TC DRESSING FORCEPS

SERRATED

51-3460	6"
51-3480	8"



TC JEWELER'S FORCEPS

DIAMOND DUST TIPS


51-4450	4"	#4
51-4451	4¼"	#5A

Sklarlite™ Forceps

Provide Lightness And Durability

SKLARLITE
THE LIGHT TOUCH

To meet industry demands for lighter pattern instruments, Sklar developed the SklarLite line of premium operating room forceps. Through a unique manufacturing technique, we are able to deliver a noticeably lighter instrument that maintains the strength and durability of a standard instrument. These extremely sensitive forceps give the surgeon maximum control and precision while comfortably withstanding the daily rigors of prolonged surgery. Sklarlite Forceps minimize fatigue and provide the exquisite control and reliability you have come to expect from all Sklar instruments. Sklarlite Forceps are covered under the same warranty that guarantees all Sklar instruments for a lifetime of dependable use.




SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ **HALSTED MOSQUITO HEMOSTATIC FORCEPS**

STRAIGHT
23-2105 5"

CURVED
23-2106 5"



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ CRILE HEMOSTATIC FORCEPS

STRAIGHT
23-2110 5½"
23-2111 7½"

CURVED
23-2113 5½"
21-2112 7½"

Sklarlite XD™


Delivers The Ultimate In Precision

SKLARLITE XD

EXTRA ◊ DELICATE

Sklar has surpassed industry standards by developing extra delicate operating room forceps that offer the surgeon one of the lightest general surgery instruments on the market today. Through a unique manufacturing technique, we are able to produce the lightest pattern forceps, affording strength and integrity equal to that of the original instrument. These extra delicate forceps are our most sensitive surgical tool, giving the surgeon the ultimate in precision and control.

Sklarlite XD Forceps are ideal for prolonged procedures, comfortably withstanding the daily rigors of frequent use. They carry the same day-in and day-out reliability that you have come to expect from all Sklar instruments and are guaranteed for a lifetime of dependable use. For the surgeon who demands the finest in OR forceps, SklarLite XD is the best choice.




SKLARLITE XD™
EXTRA ◊ DELICATE

SKLARLITE XD™ HALSTED MOSQUITO HEMOSTATIC FORCEPS

STRAIGHT
23-2505 5"

CURVED
23-2506 5"




SKLARLITE XD™
EXTRA ◊ DELICATE

SKLARLITE XD™ CRILE HEMOSTATIC FORCEPS

STRAIGHT
23-2510 5½"

CURVED
23-2513 5½"

SKLARLITE
THE LIGHT TOUCH



SKLARLITE™ KELLY HEMOSTATIC FORCEPS


STRAIGHT

23-2107 5½"

CURVED

23-2108 5½"

SKLARLITE XD™
EXTRA * DELICATE



SKLARLITE XD™ KELLY HEMOSTATIC FORCEPS


STRAIGHT

23-2507 5½"

CURVED

23-2508 5½"

SKLARLITE
THE LIGHT TOUCH



**SKLARLITE™ HARTMANN MOSQUITO
HEMOSTATIC FORCEPS**


STRAIGHT

23-2100 3½"

CURVED

23-2101 3½"

SKLARLITE XD™
EXTRA * DELICATE




**SKLARLITE XD™ HARTMANN MOSQUITO
HEMOSTATIC FORCEPS**

STRAIGHT

23-2500 3½"

CURVED

23-2501 3½"



SKLARLITE
THE LIGHT TOUCH


SKLARLITE™ ROCHESTER-PEAN FORCEPS

STRAIGHT

23-2532 6¼"
23-2534 7¼"

CURVED

23-2538 6¼"
23-2542 7¼"



SKLARLITE XD™
EXTRA * DELICATE


SKLARLITE XD™ ROCHESTER-PEAN FORCEPS

STRAIGHT

23-2552 6¼"
23-2554 7¼"

CURVED

23-2561 6¼"
23-2563 7¼"



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ ROCHESTER-OGHSNER FORCEPS

STRAIGHT 1X2 TEETH

23-2560 6¼"

CURVED 1X2 TEETH

23-2568 6¼"



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ ALLIS TISSUE FORCEPS

23-2851 6" 4x5 Teeth


23-2852 6" 5x6 Teeth



SKLARLITE XD™
EXTRA * DELICATE

SKLARLITE XD™ ALLIS TISSUE FORCEPS

23-2861	5"	4x5 Teeth
23-2862	5½"	4x5 Teeth
23-2864	6"	4x5 Teeth
23-2865	6"	5x6 Teeth



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ ADSON FORCEPS

DRESSING


23-2600	4¾"	Serrated
---------	-----	----------

TISSUE

23-2602	4¾"	1x2
---------	-----	-----

ADSON-BROWN

23-2603	4¾"	Side Grasping
---------	-----	---------------



SKLARLITE XD™
EXTRA * DELICATE

SKLARLITE XD™ ADSON FORCEPS

STRAIGHT

23-2604	4¾"	Serrated
---------	-----	----------

TISSUE

23-2606	4¾"	1x2 Teeth
---------	-----	-----------

ADSON-BROWN

23-2608	4¾"	Side Grasping
---------	-----	---------------




SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ DRESSING FORCEPS

SERRATED

23-2824	4½"
23-2826	5"
23-2828	5½"
23-2830	6"



SKLARLITE XD™
EXTRA DELICATE

SKLARLITE XD™ DRESSING FORCEPS


23-2834	4½"
23-2838	5½"
23-2836	5"
23-2840	6"



SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ TISSUE FORCEPS


1X2 TEETH	
23-2809	4½"
23-2811	5½"



SKLARLITE XD™
EXTRA DELICATE

SKLARLITE XD™ TISSUE FORCEPS

1X2 TEETH	
23-2819	5½"



HALSTED MOSQUITO HEMOSTATIC FORCEPS

STRAIGHT


17-1450 5 "

17-1650 5 " 1x2 Teeth

CURVED

17-1550 5"

17-1750 5" 1x2 Teeth




HALSTED MOSQUITO HEMOSTATIC FORCEPS

DELICATE

STRAIGHT
18-1050 5"

CURVED

18-1150 5"



CRILE HEMOSTATIC FORCEPS

STRAIGHT

17-3055 5½"

17-3062 6¼"

17-3075 7½"


CURVED

17-3155 5½"

17-3162 6¼"

17-3168 6¼" Heavy

17-3175 7½"




CRILE HEMOSTATIC FORCEPS

DELICATE

STRAIGHT
18-1952 5½"

CURVED

18-2052 5½"



KELLY HEMOSTATIC FORCEPS

STRAIGHT


17-2055 5½"

74-3260 7" Extra Heavy

CURVED

17-2155 5½"

74-3270 7" Extra Heavy



KELLY HEMOSTATIC FORCEPS


DELICATE

STRAIGHT

18-1552 5½"

CURVED

18-1652 5½"




HARTMANN MOSQUITO HEMOSTATIC FORCEPS

STRAIGHT

17-1235 3½" Serrated

CURVED

17-1335 3½" Serrated



HARTMANN MOSQUITO HEMOSTATIC FORCEPS

DELICATE

STRAIGHT

17-1230 3½"

CURVED

17-1232 3½"




ROCHESTER-PEAN FORCEPS

STRAIGHT

17-2455	5½"
17-2462	6¼"
17-2467	6¾" Heavy
17-2472	7¼"
17-2480	8"
36-2690	9"
36-2710	10¼"
36-2712	12"

CURVED

17-2255	5½"
17-2262	6¼"
17-2267	6¾" Heavy
17-2272	7¼"
17-2280	8"
36-2790	9"
36-2810	10¼"
36-2812	12"




ROCHESTER-PEAN FORCEPS

STRAIGHT, DELICATE

22-3555	5½"
---------	-----

CURVED, DELICATE

22-3655	5½"
---------	-----




ROCHESTER-OGHSNER FORCEPS

STRAIGHT, 1X2 TEETH

17-2855	5½"	(Kocher)
17-2862	6¼"	
17-2872	7¼"	
17-2880	8"	
17-2885	9"	
17-2890	10"	

CURVED, 1X2 TEETH

17-2655	5½"	(Kocher)
17-2662	6¼"	
17-2672	7¼"	
17-2680	8"	
17-2685	9"	
17-2690	10"	



ROCHESTER-OGHSNER FORCEPS


DELICATE

STRAIGHT, 1X2 TEETH

18-1752	5¼"
22-3755	5½"
17-2962	6¼"


CURVED, 1X2 TEETH

18-1852	5¼"
22-3855	5½"
17-2762	6¼"



PROVIDENCE HOSPITAL FORCEPS

STRAIGHT		CURVED	
17-1855	5½"	17-1955	5½"
18-1352	5½", Delicate	18-1452	5½", Delicate




ALLIS FORCEPS ATRAUMATIC JAWS

36-2298	6"	Straight
36-2280	8"	Straight
36-2385	10"	Straight




BOYS-ALLIS TISSUE FORCEPS

36-1860 6" 5x6 Teeth Heavy Pattern




ALLIS TISSUE FORCEPS

36-2060	6"	3x4 Teeth
36-2160	6"	4x5 Teeth
36-2260	6"	5x6 Teeth
36-2275	7½"	5x6 Teeth



ALLIS TISSUE FORCEPS

ANGLED		
36-2162	6"	4x5 Teeth
36-2262	6"	5x6 Teeth



ALLIS TISSUE FORCEPS

DELICATE		
18-2140	4"	4x5 Teeth
18-2155	5"	4x5 Teeth
22-6157	5½"	4x5 Teeth
18-2162	6"	4x5 Teeth
18-2164	6"	5x6 Teeth
36-2295	9½"	5x6 Teeth



BABCOCK FORCEPS

22-6155	5½"	6mm	Delicate
36-1662	6¼"	9mm	
36-1677	8"	10mm	
36-1692	9½"	15mm	



ADSON FORCEPS

DRESSING

50-3147	4¾"	Serrated
50-3044	4¾"	Smooth
47-2148	4¾"	Angled

TISSUE

50-3046	4¾"	Angled
50-3047	4¾"	1x2 Teeth

ADSON-BROWN

47-1447	4¾"	Side Grasping
---------	-----	---------------



ADSON FORCEPS DELICATE

DRESSING


47-2147	4¾"	Serrated
---------	-----	----------

TISSUE

47-2047	4¾"	1x2 Teeth
47-2049	4¾"	2x3 Teeth


ADSON-BROWN

47-2447	4¾"	Side Grasping
47-1450	4¾"	Angled



ADSON FORCEPS

50-3048	4¾"	Serrated, 1x2 Teeth
50-3049	4¾"	1x2 Teeth with Platform



SKLAR LITEGRIP™ ADSON FORCEPS

50-3148 Serrated 4¾"
50-3050 1x2 Teeth 4¾"
50-3052 1x2 Teeth 4¾" with Tying Platform



SKLAR LITEGRIP™ ADSON-BROWN FORCEPS

47-1448 4¾" 7x7 Teeth Side Grasping



SKLAR LITEGRIP™ BISHOP-HARMON FORCEPS

66-4034 3½" Serrated
66-4134 3½" 1x2 Teeth



DRESSING FORCEPS

SERRATED

19-1040	4"	19-1070	7"
19-1045	4½"	19-1080	8"
19-1050	5"	19-1110	10"
19-1055	5½"	19-1112	12"
19-1060	6"		



DRESSING FORCEPS

DELICATE, SERRATED

22-1161	4½"	22-1163	5½"
22-1162	5"	22-1164	6"



TISSUE FORCEPS

1X2 TEETH


19-1245	4½"	19-1280	8"
19-1250	5"	19-1310	10"
19-1255	5½"	19-1312	12"
19-1260	6"		
19-1270	7"		



TISSUE FORCEPS DELICATE

1X2 TEETH

22-1061	4½"
22-1062	5"
22-1063	5½"
22-1064	6"




TISSUE FORCEPS

2X3 TEETH

19-1445	4½"
19-1450	5"
19-1455	5½"
19-1460	6"
19-1480	8"
19-1510	10"
19-1512	12"


3X4 TEETH

19-1650	5"
19-1655	5½"
19-1660	6"
4X5 TEETH	
19-1850	5"
19-1855	5½"
19-1905	6"
19-1910	10"




SEMKIN DRESSING FORCEPS

22-4150	5"	Narrow	Serrated
22-4250	5"	Delicate	Serrated
22-4350	5"	Narrow	1x2 Teeth
22-4450	5"	Delicate	1x2 Teeth
22-4550	5"	Narrow	2x3 Teeth
22-4650	5"	Delicate	2x3 Teeth



POTTS-SMITH FORCEPS


52-3270	7"	Serrated	52-3410	10"	1x2 Teeth
52-3370	7"	1x2 Teeth	52-3510	10"	1x2 Teeth, Curved
52-3282	8¼"	Serrated	52-3312	12"	Serrated
52-3382	8¼"	1x2 Teeth	52-3412	12"	1x2
52-3310	10"	Serrated			



DEBAKEY ATRAUMATIC TISSUE FORCEPS

1.5MM JAWS


52-5062	6¼"
52-5077	7¾"
52-5095	9½"
52-5098	12"



DEBAKEY ATRAUMATIC TISSUE FORCEPS

2MM JAWS


52-5162	6¼"
52-5177	7¾"
52-5195	9½"
52-5198	12"



DEBAKEY ATRAUMATIC TISSUE FORCEPS

3.5MM JAWS


52-5262	6¼"
52-5277	7¾"
52-5295	9½"
52-5297	12"



EYE DRESSING FORCEPS

4", SERRATED

66-2734	Straight	.5mm Delicate
66-2736	Half Curve	.5mm Delicate
66-2738	Full Curve	.5mm Delicate
66-2740	Straight	.8mm
66-2840	Straight	.8mm Non-Magnetic
66-2940	Half Curve	.8mm
66-3140	Full Curve	.8mm




EYE TISSUE FORCEPS

1X2 TEETH, 4"

66-3340	Straight	.8mm
66-3240	Straight	.8mm, Non-Magnetic
66-3540	Half Curve	.8mm
66-3740	Full Curve	.8mm

1X2 TEETH, EXTRA DELICATE, 4"

66-3341	Straight	.5mm
66-3541	Half Curve	.5mm
66-3741	Full Curve	.5mm



GRAEFE IRIS FORCEPS

STRAIGHT

66-4526	2¾"	Serrated
66-4525	2¾"	1x2 Teeth


CURVED

66-4626	2¾"	Serrated
66-4625	2¾"	1x2 Teeth



BISHOP-HARMON FORCEPS

66-4032	3½"	Serrated
66-4132	3½"	1x2 Teeth
66-4230	3½"	1x2 Teeth Micro
66-4232	3½"	1x2 Teeth Delicate
66-4133	3½"	1x2 Teeth Extra Delicate



GRAEFE FIXATION FORCEPS

66-1940 4" Without Lock
66-2140 4" With Lock



HARMON FORCEP

66-1930 3"



JEWELERS FORCEPS

66-7436 4½" With Extra Fine Tips,
Style 1



JEWELERS FORCEPS

66-7438 4½" With Fine Tips, Style 1



JEWELERS FORCEPS

66-7440 4½" With Fine Tips, Style 2




**JEWELERS STYLE
FORCEPS**

66-7447 4½" #3 Straight
66-7442 4½" #3C Straight, Delicate
66-7443 4¾" #3 Straight, Micro-Fine



JEWELERS STYLE FORCEPS

DELICATE, STRAIGHT
66-7542 4½" #4



JEWELERS STYLE FORCEPS

STRAIGHT
66-7642 4½" #5



JEWELERS STYLE FORCEPS

ANGLED
66-7645 4½" #5




JEWELERS STYLE FORCEPS

CURVED
66-7644 4½" #7



FINE POINT SPLINTER FORCEPS

STRAIGHT
19-3035 3½"
19-3040 4"
19-3045 4½"



CARMALT SPLINTER FORCEPS

19-3240 4¼" Straight
19-3340 4¼" Curved




WALTER SPLINTER FORCEPS

19-3250 4¼" Straight
19-3252 4¼" Curved



FEILCHENFELD SPLINTER FORCEPS


STRAIGHT
19-2930 3"
19-2945 4½"



PLAIN SPLINTER FORCEPS

WITH MAGNIFIER


19-3030 3" Non-Autoclavable



**VANGUARD SPLINTER
LIBERATOR**


TROCAR POINT

19-3350 3¼" Handle, 4½" Overall



**SKLAR SPLINTER
REMOVAL KIT**


19-3360 Containing one 19-3240 and
one 19-3350 in a Convenient
Storage Case.



VIRTUS SPLINTER FORCEPS

SERRATED


19-3160 6" Straight
19-3162 6" Angled



ARCHER SPLINTER FORCEPS

SERRATED


19-3155 5½" Straight



PEET SPLINTER FORCEPS

DOUBLE EDGE

19-3147 4¾" Angled



WILLIAMS (STIEGLITZ) SPLINTER FORCEPS

19-3056	4¾"	Right Angle
19-3055	5½"	45° Angle
19-3161	5½"	Straight



ZIEGLER CILIA FORCEPS

SMOOTH		
66-1535	3½"	2mm




BERGH CILIA FORCEPS

66-1635	3½"	5mm
---------	-----	-----



SWISS CILIA FORCEPS

DIAGONAL JAWS		
66-1740	3¾"	3mm



DOUGLASS CILIA FORCEPS

THIN TIPS		
66-1830	3½"	2.5mm



LITTAUER CILIA FORCEPS

66-1735	3½"	4mm
---------	-----	-----



BARRAQUER CILIA FORCEPS

SMOOTH

66-1742 4¼" 5mm



BEER CILIA FORCEPS


SMOOTH

66-1770 3½" 4mm




SHAAF CILIA FORCEPS

66-6034 3½"



DESMARRES CHALAZION FORCEPS

66-1038	3½"	Small
66-1036	3½"	Medium
66-1035	3½"	Large




BALLENGER SPONGE FORCEPS

STRAIGHT

74-1370	7"	Smooth
74-1470	7"	Serrated

CURVED

74-1570	7"	Smooth
74-1670	7"	Serrated




FOERSTER SPONGE FORCEPS

STRAIGHT

87-2095	9½"	Smooth
87-2195	9½"	Serrated

CURVED


87-2295	9½"	Smooth
87-2395	9½"	Serrated



FLETCHER SPONGE FORCEPS


SERRATED

93-1222	9½"	Straight
93-1226	9½"	Curved




BACKHAUS TOWEL CLAMPS

47-2935	3½"
47-2945	4½"
47-2955	5¼"
47-2960	6"



JONES CROSS ACTION TOWEL CLAMPS

47-2925	2¼"	Delicate
47-2927	3"	
47-2930	3½"	




STERILIZATION FORCEPS

STRAIGHT

10-1080	8"
10-1111	11"

CURVED


10-1081	8"
10-1211	11"



KERN BONE HOLDING FORCEPS

WITHOUT RATCHET


40-2480	5¾"	Extra Small
40-2485	8½"	



KERN BONE HOLDING FORCEPS


With Ratchet

40-2580	5¾"	Extra Small
40-2585	8½"	




BONE AND CARTILAGE CLAMP

40-2275 8" Curved



MARTIN CARTILAGE CLAMP


40-2273 7½"



**DINGMAN SMALL BONE AND
CARTILAGE CLAMP**


SERRATED JAWS

40-2660 7½" 4mm 2x2 Teeth




SESAMOID CLAMP

40-2655 6½"



LOCKE PHALANGEAL FORCEPS

40-2650 6"



LEWIN BONE FORCEPS


40-2670 7" Serrated



LOWMAN CLAMPS

1X2 PRONG


40-2921	4¾"x¾" Jaw
40-2922	7¼"x1½" Jaw
40-2923	8"x1½" Jaw



LAMBERT-LOWMAN BONE CLAMP

2X2 PRONG

40-2838	7" 1" Jaw
---------	-----------



LOWMAN-HOGLUND BONE CLAMP

2X2 PRONG, SWIVEL HEAD

40-2828	7½" 1" Jaw
---------	------------



BOZEMAN FORCEPS


ONE LARGE RING

90-7710	10¼" S-Curve
---------	--------------




BOZEMAN FORCEPS

90-7610	10½"	Straight
80-1410	10¼"	S-Curve



GRASPING FORCEPS


10-1013	6¼"
---------	-----



**CARROLL TENDON PULLING
FORCEP**


ALLIGATOR TYPE

47-3520	4¼"	Curved
---------	-----	--------



**MESH BUR AND DIAMOND
HOLDERS**


49-3155	1⅜"
49-3150	2⅛"



SKLAR EDGE™ TC OPERATING SCISSORS

STRAIGHT

16-1055	5½"	Sharp/Sharp
16-1056	5½"	Sharp/Blunt
16-1057	5½"	Blunt/Blunt



SKLAR EDGE™ TC OPERATING SCISSORS

CURVED

16-1155	5½"	Sharp/Sharp
16-1156	5½"	Sharp/Blunt
16-1157	5½"	Blunt/Blunt



SKLAR EDGE™ TC IRIS SCISSORS

SHARP/SHARP


16-2405	4½"	Straight
16-2505	4½"	Curved



SKLAR EDGE™ TC MAYO DISSECTING SCISSORS

BLUNT/BLUNT


16-1505	5½"	Straight
16-1605	5½"	Curved
16-1510	6¾"	Straight
16-1610	6¾"	Curved
16-1515	9"	Straight
16-1615	9"	Curved



SKLAR EDGE™ TC MAYO-STILLE DISSECTING SCISSORS

BLUNT/BLUNT

16-3552	5½"	Straight
16-3562	5½"	Curved
16-3557	6¾"	Straight
16-3567	6¾"	Curved



SKLAR EDGE™ TC METZENBAUM SCISSORS (LAHEY)

BLUNT/BLUNT

16-1810	5¾"	Straight	
16-1910	5¾"	Curved	
16-1805	5¾"	Straight	Delicate
16-1905	5¾"	Curved	Delicate

SHARP/SHARP

16-1815	5¾"	Straight
16-1915	5¾"	Curved



**SKLAR EDGE™ TC
METZENBAUM SCISSORS**

STRAIGHT, BLUNT/BLUNT

16-2010	7"	
16-2005	7"	Delicate
16-2210	9"	(Nelson)
16-2205	9"	Delicate (Nelson)

STRAIGHT, SHARP/SHARP

16-2015	7"	Delicate
---------	----	----------



**SKLAR EDGE™ TC
METZENBAUM SCISSORS**

CURVED, BLUNT/BLUNT

16-2110	7"	
16-2105	7"	Delicate
16-2215	8"	
16-2310	9"	
16-2305	9"	Delicate (Nelson)
16-2211	10"	Delicate (Nelson)
16-2311	11"	(Nelson)
16-2312	12"	(Nelson)


CURVED, SHARP/SHARP

16-2115	7"	Delicate
---------	----	----------



**SKLAR EDGE™ TC STRABISMUS
SCISSORS**


16-2605	4½"	Straight
16-2705	4½"	Curved



**SKLAR EDGE™ TC
GOLDMAN-FOX
SCISSORS**


ONE SERRATED BLADE

16-3614	5"	Straight	Blunt
16-3615	5"	Curved	Blunt
16-2610	5"	Curved	Sharp




**SKLAR EDGE™ TC
KAYE RHYTIDECTOMY
& DISSECTING
SCISSORS**

16-3621	4½"	Curved
16-3620	5¾"	Curved



**SKLAR EDGE™ TC WIRE
CUTTING SCISSORS**

16-2800	4¾"	Serrated	Angled
---------	-----	----------	--------




**SKLAR EDGE™ TC
STEVENS TENOTOMY
SCISSORS**

BLUNT/BLUNT


16-2730	4½"	Straight
16-2731	4½"	Curved

SUPER-CUT SKLAR EDGE™ TUNGSTEN CARBIDE SKLARHONE™ scissors combine all the advantages of durable tungsten carbide with the precision performance of SUPER-CUT SklarHone™ blades, featuring the razor sharp cutting edge of a knife on one blade and micro-serrations on the other. Easily identified with one gold handle and one black handle.




SUPER-CUT TC SKLARHONE™ IRIS SCISSORS

16-3545	4½"	Straight
16-3546	4½"	Curved



SUPER-CUT TC SKLARHONE™ MAYO SCISSORS

16-3530	5½"	Straight
16-3531	5½"	Curved
16-3532	6¾"	Straight
16-3533	6¾"	Curved



SUPER-CUT TC SKLARHONE™ MAYO-STILLE SCISSORS

16-3572	6"	Straight
16-3582	6"	Curved
16-3577	6¾"	Straight
16-3587	6¾"	Curved



SUPER-CUT TC SKLARHONE™ METZENBAUM SCISSORS

STRAIGHT, BLUNT/BLUNT TIPS


16-3503	5¾"	(Lahey)
16-3513	5¾"	Delicate (Lahey)
16-3504	7"	
16-3514	7"	Delicate
16-3505	8"	
16-3515	8"	Delicate



SUPER-CUT TC SKLARHONE™ METZENBAUM SCISSORS

CURVED, BLUNT/BLUNT TIPS


16-3506	5¾"	(Lahey)
16-3526	5¾"	Delicate (Lahey)
16-3507	7"	
16-3527	7"	Delicate
16-3508	8"	
16-3528	8"	Delicate
16-3509	9"	(Nelson)
16-3510	10"	(Nelson)
16-3511	11"	(Nelson)



**SUPER-CUT TC SKLARHONE™
STEVENS TENOTOMY SCISSORS**


POINTED, SHARP/SHARP

16-3540 4½" Straight
16-3541 4½" Curved



SUPER-CUT TC SKLARHONE™ KILNER SCISSORS

16-3619 5" Curved



**SKLARHONE™
OPERATING SCISSORS**

CURVED, SMOOTH


15-3349 5½" Sharp/Blunt
15-3348 5½" Sharp/Sharp

**SUPER-CUT
SKLARHONE™**

**SUPER-CUT SKLARHONE™
OPERATING SCISSORS**

CURVED, SERRATED

15-3449 5½" Sharp/Blunt
15-3448 5½" Sharp/Sharp



**SKLARHONE™
IRIS SCISSORS**

SMOOTH

15-3352 4½" Straight
15-3357 4½" Curved

**SUPER-CUT
SKLARHONE™**

**SUPER-CUT SKLARHONE™
IRIS SCISSORS**

SERRATED

15-3452 4½" Straight
15-3457 4½" Curved

SKLARHONE™ MAYO SCISSORS

SMOOTH

15-3325	5½"	Straight
15-3330	5½"	Curved
15-3326	6¾"	Straight
15-3331	6¾"	Curved

SUPER-CUT SKLARHONE™

SUPER-CUT SKLARHONE™ MAYO SCISSORS

SERRATED

15-3429	5½"	Straight
15-3430	5½"	Curved
15-3426	6¾"	Straight
15-3427	6¾"	Curved
15-3434	9"	Straight
15-3435	9"	Curved



SKLARHONE™ METZENBAUM SCISSORS

SMOOTH


15-3305	5¾"	Straight	(Lahey)
15-3306	5¾"	Curved	(Lahey)
15-3507	7"	Straight	
15-3310	7"	Curved	
15-3508	8"	Straight	
15-3311	8"	Curved	
15-3509	9"	Straight	(Nelson)
15-3315	9"	Curved	(Nelson)

SUPER-CUT SKLARHONE™

SUPER-CUT SKLARHONE™ METZENBAUM SCISSORS

SERRATED

15-3405	5¾"	Straight	(Lahey)
15-3406	5¾"	Curved	(Lahey)
15-3407	7"	Straight	
15-3410	7"	Curved	
15-3408	8"	Straight	
15-3411	8"	Curved	
15-3409	9"	Straight	(Nelson)
15-3412	9"	Curved	(Nelson)



SKLARHONE™ SKLARLITE™ METZENBAUM SCISSORS

SMOOTH

15-9107	7"	Straight	
15-9117	7"	Curved	
15-9108	8"	Straight	
15-9118	8"	Curved	
15-9109	9"	Straight	(Nelson)
15-9119	9"	Curved	(Nelson)


SKLARLITE THE LIGHT TOUCH

SUPER-CUT SKLARHONE™

SUPER-CUT SKLARHONE™ SKLARLITE™ METZENBAUM SCISSORS

SERRATED

15-3307	7"	Straight	
15-3317	7"	Curved	
15-3308	8"	Straight	
15-3318	8"	Curved	
15-3309	9"	Straight	(Nelson)
15-3319	9"	Curved	(Nelson)




SUPER-CUT SKLARHONE™

SUPER-CUT SKLARHONE™ GOLDMAN-FOX SCISSOR

SERRATED


15-3340	5"	Curved
---------	----	--------



SUPER-CUT
SKLARHONE™

**SUPER-CUT SKLARHONE™
GRADLE SCISSOR**


SERRATED
15-3444 3¾" Slight Curve



SUPER-CUT
SKLARHONE™

**SUPER-CUT SKLARHONE™
STEVENS TENOTOMY
SCISSORS**


SERRATED
15-1242 4½" Straight Blunt
15-1243 4½" Straight Sharp
15-1342 4½" Curved Blunt
15-1343 4½" Curved Sharp



SUPER-CUT
SKLARHONE™

**SUPER-CUT
SKLARHONE™
LISTER SCISSORS**

SERRATED
15-3355 5½"
15-3356 7¼"
15-3358 8"






SKLARCUT™ OPERATING SCISSORS

15-3505	5½"	Straight	Sharp/Blunt
15-3500	5½"	Straight	Sharp/Sharp
15-3515	5½"	Curved	Sharp/Blunt
15-3510	5½"	Curved	Sharp/Sharp



SKLARLITE™ SKLARCUT™ OPERATING SCISSORS

23-1262	5½"	Straight	Sharp/Blunt
23-1260	5½"	Straight	Sharp/Sharp
23-1266	5½"	Curved	Sharp/Blunt
23-1264	5½"	Curved	Sharp/Sharp




SKLARCUT™ IRIS SCISSORS

15-3585	4½"	Straight
15-3590	4½"	Curved




SKLARLITE™ SKLARCUT™ IRIS SCISSORS

23-1268	4½"	Straight
23-1270	4½"	Curved



SKLARCUT™ STEVENS TENOTOMY SCISSORS

BLUNT/BLUNT		
15-3372	4½"	Straight
15-3370	4½"	Curved




SKLARCUT™ MAYO SCISSORS

15-3550	5½"	Straight
15-3560	5½"	Curved
15-3555	6¾"	Straight
15-3565	6¾"	Curved
15-3566	9"	Straight
15-3568	9"	Curved




**SKLARLITE™ SKLARCUT™
MAYO SCISSORS**

23-1274 6¾" Straight
23-1278 6¾" Curved



**SKLARCUT™ METZENBAUM
SCISSORS (LAHEY)**

15-3569 4½" Straight
15-3571 4½" Curved
15-3567 5¾" Straight
15-3570 5¾" Curved



**SKLARLITE™ SKLARCUT™
METZENBAUM SCISSORS**


23-1280 5¾" Curved (Lahey)
23-1267 7" Curved



**SKLARCUT™
METZENBAUM
SCISSORS**


BLUNT/BLUNT

15-3575 7" Curved
15-3586 7" Straight
15-3574 9" Curved




**SKLARCUT™
GOLDMAN-FOX
SCISSORS**

15-3536 5" Straight
Sharp/Sharp



**SKLARCUT™ WAGNER
SCISSORS**


15-3542 4½" Straight
Sharp/Blunt



**SKLARLITE™
SKLARCUT™ WAGNER
SCISSORS**


23-1254 4½" Straight
Sharp/Blunt

SKLARCUT
SERRATED BLADE




SKLARCUT™ GRADLE SCISSORS

15-3244 3¾" Curved Blunt, Serrated




SKLARCUT™ KAYE FACE LIFT SCISSORS

15-3579 6" Curved Blunt/Blunt
15-3580 7½" Curved Blunt/Blunt



SKLARCUT™ FOMON LATERAL SCISSORS

15-3548 5" Curved Blunt/Blunt



SKLARCUT™ COTTLE DORSAL SCISSORS

15-3547 6½" Angled Blunt/Blunt



SKLARCUT™ CASTANARES FACE LIFT SCISSOR

15-3549 6¼" Straight Blunt/Blunt

SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™ OPERATING
SCISSORS**

STRAIGHT


23-1104	4½"	Sharp/Blunt
23-1100	4½"	Sharp/Sharp
23-1105	5"	Sharp/Blunt
23-1101	5"	Sharp/Sharp
23-1106	5½"	Sharp/Blunt
23-1102	5½"	Sharp/Sharp

SKLARLITE XD
EXTRA • DELICATE

**SKLARLITE XD™ OPERATING
SCISSORS**

STRAIGHT

23-1116	4½"	Sharp/Blunt
23-1112	4½"	Sharp/Sharp
23-1117	5"	Sharp/Blunt
23-1113	5"	Sharp/Sharp
23-1118	5½"	Sharp/Blunt
23-1114	5½"	Sharp/Sharp



SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™ OPERATING
SCISSORS**

CURVED


23-1128	4½"	Sharp/Blunt
23-1124	4½"	Sharp/Sharp
23-1130	5½"	Sharp/Blunt
23-1126	5½"	Sharp/Sharp

SKLARLITE XD
EXTRA • DELICATE

**SKLARLITE XD™ OPERATING
SCISSORS**

CURVED

23-1140	4½"	Sharp/Blunt
23-1136	4½"	Sharp/Sharp
23-1142	5½"	Sharp/Blunt
23-1138	5½"	Sharp/Sharp



SKLARLITE
THE LIGHT TOUCH


**SKLARLITE™ MAYO
DISSECTING SCISSORS**

23-1164	5½"	Straight
23-1166	5½"	Curved
23-1165	6¾"	Straight
23-1167	6¾"	Curved

SKLARLITE XD
EXTRA • DELICATE

**SKLARLITE XD™ MAYO
DISSECTING SCISSORS**

23-1168	5½"	Straight
23-1170	5½"	Curved
23-1169	6¾"	Straight
23-1171	6¾"	Curved



SKLARLITE
THE LIGHT TOUCH


**SKLARLITE™ METZENBAUM
SCISSORS (LAHEY)**

23-1172	4½"	Straight
23-1174	4½"	Curved
23-1173	5¾"	Straight
23-1175	5¾"	Curved

SKLARLITE XD
EXTRA • DELICATE

**SKLARLITE XD™ METZENBAUM
SCISSORS (LAHEY)**

23-1176	4½"	Straight
23-1178	4½"	Curved
23-1177	5¾"	Straight
23-1179	5¾"	Curved



SKLARLITE™ THE LIGHT TOUCH

SKLARLITE™ METZENBAUM SCISSORS

23-1180	7"	Straight
23-1182	7"	Curved



SKLARLITE XD™ EXTRA * DELICATE

SKLARLITE XD™ METZENBAUM SCISSORS


23-1184	7"	Straight
23-1186	7"	Curved

SKLARLITE™ THE LIGHT TOUCH

SKLARLITE™ IRIS SCISSORS

SHARP/SHARP

23-1148	3½"	Straight
23-1149	4½"	Straight
23-1150	3½"	Curved
23-1151	4½"	Curved



SKLARLITE XD™ EXTRA * DELICATE

SKLARLITE XD™ IRIS SCISSORS


SHARP/SHARP

23-1152	3½"	Straight
23-1153	4½"	Straight
23-1154	3½"	Curved
23-1155	4½"	Curved

SKLARLITE™ THE LIGHT TOUCH

SKLARLITE™ STEVENS SCISSORS

23-1156	4½"	Straight	Blunt/Blunt
23-1157	4½"	Curved	Blunt/Blunt



SKLARLITE XD™ EXTRA * DELICATE

SKLARLITE XD™ STEVENS SCISSORS


BLUNT/BLUNT

23-1160	4½"	Straight	Blunt/Blunt
23-1161	4½"	Curved	Blunt/Blunt

SKLARLITE™ THE LIGHT TOUCH

SKLARLITE™ STITCH SCISSORS

23-1188	3½"	Spencer
23-1189	4½"	Spencer
23-1190	5½"	Littauer



SKLARLITE XD™ EXTRA * DELICATE


SKLARLITE XD™ STITCH SCISSORS

23-1193	3½"	Spencer
23-1194	4½"	Spencer
23-1195	5½"	Littauer

SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ SHORTBENT
SUTURE SCISSORS

23-1191 3½"



SKLARLITE XD
EXTRA • DELICATE

SKLARLITE XD™ SHORTBENT
SUTURE SCISSOR

23-1196 3½"

SKLARLITE
THE LIGHT TOUCH

SKLARLITE™ LISTER
BANDAGE SCISSORS

23-1199 5½"

23-1200 7¼"



SKLARLITE XD
EXTRA • DELICATE

SKLARLITE XD™ LISTER
BANDAGE SCISSORS

23-1202 5½"

23-1203 7¼"

SKLAR DISINFECTANT™ SURFACE WIPES

Bactericidal • Virucidal • Fungicidal • Tuberculocidal

Ready-To-Use / Contains Biodegradable Detergent / Non-Bleach Formula

Kills VRE and MRSA


2-Minute Kill Time!

Item# 10-1616 Case of 6 Canisters

Item# 10-1617 Case of 12 Canisters

100 Towelettes per Canister





OPERATING SCISSORS

STRAIGHT


SHARP/SHARP		SHARP/BLUNT		BLUNT/BLUNT	
13-1045	4½"	14-1045	4½"	15-1045	4½"
13-1050	5"	14-1050	5"	15-1050	5"
13-1055	5½"	14-1055	5½"	15-1055	5½"
13-1060	6"	14-1060	6"	15-1060	6"
13-1065	6½"	14-1065	6½"	15-1065	6½"
13-1072	7¼"	14-1072	7¼"		



OPERATING SCISSORS

DELICATE, STRAIGHT


SHARP/SHARP		SHARP/BLUNT		BLUNT/BLUNT	
22-1145	4½"	22-1245	4½"	22-1345	4½"
22-1150	5"	22-1250	5"	22-1350	5"
22-1155	5½"	22-1255	5½"	22-1355	5½"



OPERATING SCISSORS

CURVED


SHARP/SHARP		SHARP/BLUNT		BLUNT/BLUNT	
13-2045	4½"	14-2045	4½"	15-2045	4½"
13-2050	5"	14-2050	5"	15-2050	5"
13-2055	5½"	14-2055	5½"	15-2055	5½"
13-2060	6"	14-2060	6"	15-2060	6"
13-2065	6½"	14-2065	6½"	15-2065	6½"
13-2072	7¼"	14-2072	7¼"		



OPERATING SCISSORS

DELICATE, CURVED

SHARP/SHARP		SHARP/BLUNT		BLUNT/BLUNT	
22-1445	4½"	22-1545	4½"	22-1645	4½"
22-1450	5"	22-1550	5"	22-1650	5"
22-1455	5½"	22-1555	5½"	22-1655	5½"



DEAVER SCISSORS


DELICATE

STRAIGHT

13-1155	5½"	Sharp/Sharp
14-1155	5½"	Sharp/Blunt
15-1155	5½"	Blunt/Blunt

CURVED

13-2155	5½"	Sharp/Sharp
14-2155	5½"	Sharp/Blunt



WAGNER SCISSORS


STRAIGHT

SMOOTH

47-1047	4½"	Sharp/Sharp
47-1048	4½"	Sharp/Blunt
47-1070	4½"	Blunt/Blunt

SERRATED

47-1049	4½"	Sharp/Sharp
47-1060	4½"	Sharp/Blunt
47-1071	4½"	Blunt/Blunt



WAGNER SCISSORS


CURVED

SMOOTH

47-1058	4½"	Sharp/Sharp
47-1068	4½"	Sharp/Blunt
47-1075	4½"	Blunt/Blunt

SERRATED

47-1059	4½"	Sharp/Sharp
47-1069	4½"	Sharp/Blunt
47-1072	4½"	Blunt/Blunt




MAYO DISSECTING SCISSORS

STRAIGHT

15-1555	5½"
15-1567	6¾"
90-1590	9"
55-8711	11" Harrington

CURVED

15-2555	5½"
15-2567	6¾"
90-1690	9"
55-8611	11" Harrington



MAYO DISSECTING SCISSORS


DELICATE

STRAIGHT

22-1556	5½"
22-1567	6¾"

CURVED

22-2555	5½"
22-2567	6¾"




METZENBAUM (LAHEY) SCISSORS

STRAIGHT

22-2550	4½"
75-5460	5¾"

CURVED

22-2545	4½"
47-1360	5¾"



METZENBAUM (LAHEY) SCISSORS


DELICATE

STRAIGHT

22-2450	4½"
22-2460	5¾"

CURVED

22-2650	4½"
22-1057	5¾"




METZENBAUM DISSECTING SCISSORS

STRAIGHT

75-5570	7"
55-9180	8" (Nelson)
55-9190	9" (Nelson)
55-9310	10" (Nelson)
55-9311	11" (Nelson)
55-9312	12" (Nelson)
55-9314	14" (Nelson)

CURVED

75-5670	7"
55-9080	8" (Nelson)
55-9090	9" (Nelson)
55-9210	10" (Nelson)
55-9411	11" (Nelson)
55-9612	12" (Nelson)
55-9614	14" (Nelson)



METZENBAUM DISSECTING SCISSORS


DELICATE

STRAIGHT

22-1071	7"	
55-9184	8"	(Nelson)
55-9194	9"	(Nelson)

CURVED

22-1070	7"	
55-9084	8"	(Nelson)
55-9094	9"	(Nelson)
55-9414	11"	(Nelson)




LEXER SCISSORS

15-2426	6¼"	Straight
15-2436	6¼"	Curved
15-2428	8¼"	Straight
15-2438	8¼"	Curved

BABY LEXER

22-1040	4"	Straight
22-1042	4"	Curved



MAYO-STILLE DISSECTING SCISSORS


ROUNDED BLADES

STRAIGHT

15-1655	5½"
15-1667	6¾"

CURVED


15-2655	5½"
15-2667	6¾"



KELLY-ADSON SCISSORS

STRAIGHT


50-4562	6¼"	Smooth
50-4563	6¼"	Serrated



KELLY-ADSON SCISSORS

CURVED


50-4662 6¼" Smooth



KELLY-ADSON SCISSORS

ANGULAR

50-4462 6¼" Smooth



IRIS SCISSORS

64-2142	4½"	Angled
64-2143	4½"	Curved, Sideways
64-2042	4¼"	Angled, Probe Point



IRIS SCISSORS

STRAIGHT

47-1135	3½"	S/S
47-1140	3½"	B/B
47-1143	4"	S/S
47-1441	4"	S/S, Heavy
47-1145	4½"	S/S
47-1144	4½"	B/B

CURVED

47-1235	3½"	S/S
47-1240	3½"	B/B
47-1243	4"	S/S
47-1442	4"	S/S, Heavy
47-1245	4½"	S/S
47-1244	4½"	B/B



IRIS SCISSORS DELICATE


SHARP/SHARP

STRAIGHT

22-1135	3½"
47-1440	4"
47-1246	4½"

CURVED

22-1235	3½"
47-1340	4"
47-1248	4½"




STRABISMUS SCISSORS

STRAIGHT

64-1442	4½"
---------	-----

CURVED

64-1542	4½"
---------	-----



STRABISMUS SCISSORS


DELICATE

STRAIGHT

22-1442	4½"
---------	-----

CURVED

22-1542	4½"
---------	-----



STEVENS TENOTOMY SCISSORS

STRAIGHT

64-1243	4½"	S/S
64-1242	4½"	B/B

STRAIGHT, DELICATE


22-1240	4½"	S/S
22-1242	4½"	B/B

CURVED

64-1343	4½"	S/S
64-1342	4½"	B/B

CURVED, DELICATE


22-1241	4½"	S/S
22-1342	4½"	B/B



KAYE UNDERMINING SCISSORS

STRAIGHT

47-1084 4½"




RIBBON STRABISMUS SCISSORS

STRAIGHT

64-1240 4"

CURVED

64-1340 4"




RIBBON STEVENS TENOTOMY SCISSORS

STRAIGHT

64-1241 4"

CURVED

64-1341 4"



BONN MINI IRIS SCISSORS

STRAIGHT

64-2035 3½"

CURVED


64-2036 3½"



REYNOLDS SCISSORS


NARROW TIP

22-2758	6"	Straight
22-2759	6"	Straight, Serrated
22-2760	6"	Curved
22-2762	6"	Curved, Serrated
22-2772	7"	Straight
22-2768	7"	Straight, Serrated
22-2770	7"	Curved
22-2722	7"	Curved, Serrated



**GOLDMAN-FOX
SCISSORS**


CURVED
22-2520 5" Serrated



JOSEPH SCISSORS

DOUBLE EDGE

STRAIGHT **CURVED**
41-1058 5¾" 41-1060 5¾"




JOSEPH SCISSORS

STRAIGHT

41-1050 5¾"


CURVED

41-1054 5¾"




**KAYE FACELIFT
SCISSORS**

47-1086 6"
47-1087 7½"




FOMON SCISSORS

70-4352 5¼" Curved
41-1076 5½" Curved,
Saber Back
70-4350 5½" Angled




**CASTANARES
FACE LIFT SCISSORS**

41-1077 6¼" Straight



**COTTLE DORSAL
SCISSORS**

HEAVY ROUNDED BLADES
41-1169 6" Angular




**COTTLE NASAL
SCISSORS**

STRONG SHORT BLADES
41-1168 4½" Curved




**RAGNELL DISSECTING
SCISSORS (FLAT TIPS)**

CURVED	STRAIGHT
22-2652 5"	22-2653 5"
22-2657 7"	




**KILNER (RAGNELL)
DISSECTING
SCISSORS**

CURVED
47-1042 4¾" Blunt/Blunt
Smooth




**KILNER DISSECTING
SCISSORS**

CURVED FLAT TIPS
47-1043 6" Blunt/Blunt




WILMER IRIS SCISSORS

ANGLED ON FLAT
64-1735 4"




SISTRUNK SCISSORS

15-2445 5¾"	Straight
15-2455 5¾"	Curved



LA GRANGE SCISSORS

S-SHAPED
22-2745 4½" Sharp/Sharp




**SKLAR NAIL
SCISSORS**

RING HANDLE
97-1235 3½" Straight
97-1245 4½" Curved




**BIRO DERMAL
NAEVUS SCISSORS**

DELICATE
64-1740 4"



**BROPHY SULLIVAN
SCISSORS**

22-2755 5¾" Serrated Curved
22-2855 5¾" Serrated Straight



STITCH SCISSORS

22-2835 3½" Spencer
22-2845 4½" Spencer
15-1355 5½" Littauer


DELICATE

22-2834 3½" Spencer
22-2844 4½" Spencer
22-1354 5½" Littauer



SUTURE SCISSORS

22-2935 3½" Shortbent
22-2934 3½" Shortbent, Delicate
22-2947 4¾" Northbent




**SKLAR STITCH
SCISSORS**

15-3150 4½" Angled
22-2950 5" Curved




HILL STITCH SCISSORS

22-2955 4¾" Curved




**O'BRIEN STITCH
SCISSORS**

22-2837 3¾" Angled




**SUTURE WIRE
CUTTING SCISSORS**

24-2342 4¾" Angled,
Serrated




**LONG STITCH
SCISSORS**

15-1358 6¼"



BUCK SUTURE SCISSORS

15-1359 7"




GRADLE SCISSORS

64-3280 3¾" Slight Curve, Sharp
22-2957 3¾" Ribbon Style



CROWN SCISSORS


24-2345 4" Straight, Smooth
24-2340 4" Straight, Serrated
24-2346 4" Curved, Smooth
24-2341 4" Curved, Serrated



MOLESKIN SCISSORS

ONE BLADE SERRATED, SHARP/BLUNT


97-1350 7½" Straight
97-1351 7½" Curved



INGROWN NAIL SPLITTING SCISSORS


ONE POINTED BLADE

97-0555 6" One Blade Serrated



LISTER BANDAGE SCISSORS


11-1035 3½"
11-1045 4½"
11-1055 5½"
11-1055LFT 5½" Left Handed
11-1072 7¼"
11-1072LFT 7¼" Left Handed
11-1080 8"



LISTER BANDAGE SCISSORS

DELICATE


22-1054 5½"
22-1072 7¼"



LISTER BANDAGE SCISSORS

LARGE RING


11-1060 8"



SKLAR HI LEVEL BANDAGE SCISSORS


SERRATED

11-1252 5½"




NURSES BANDAGE SCISSORS

11-1495 5½" Pink
11-1496 5½" Blue




KNOWLES BANDAGE SCISSORS

11-1255 5½"



USA PATTERN GAUZE SCISSORS


11-1180 8"



UTILITY SCISSORS

SERRATED


11-1278	6"	Blue
11-1280	7½"	Blue
11-1281	7½"	Black
11-1282	7½"	Green
11-1283	7½"	Red
11-1284	7½"	White
11-1285	7½"	Purple
11-1286	7½"	Orange
11-1287	7½"	Yellow



BRUNS PLASTER SHEARS

HEAVY PATTERN

11-1380	8"
11-1390	9"



HERCULES UTILITY SCISSORS


11-1297	7½"	Serrated
---------	-----	----------



SEUTIN PLASTER SHEARS

HEAVY PLIER HANDLE


11-1490	9"	Serrated
---------	----	----------



TC DERF NEEDLE HOLDER

SERRATED

21-8004 4¾"




TC WEBSTER NEEDLE HOLDERS

21-2240 4¼" Smooth

21-8002 5¼" Smooth


21-8080 5¼" Serrated



TC HALSEY NEEDLE HOLDERS

21-8000 5¼" Smooth


21-8001 5¼" Serrated



TC NEIVERT NEEDLE HOLDER

SERRATED


21-8086 5" One Offset Ring



TC BAUMGARTNER NEEDLE HOLDER

SERRATED

21-8005 5¼"




TC PAR NEEDLE HOLDERS

DELICATE, SERRATED

21-7012 4½"

21-7013 5"



TC CRILE-WOOD NEEDLE HOLDERS

21-8003 6" Smooth

21-8006 6" Serrated


21-8007 7" Serrated

21-8008 8" Serrated

21-8009 9" Serrated

21-8010 10½" Serrated

21-8011 12" Serrated



TC MAYO-HEGAR NEEDLE HOLDERS

SERRATED

21-8030 5½"

21-8031 6"

21-8032 7"

21-8033 8"

21-8051 9½"

21-8034 10½"


21-8035 12"



TC RYDER NEEDLE HOLDERS

MICRO SERRATED


21-8012	5"
21-8013	6"
21-8014	7"
21-8016	9"



TC RYDER NEEDLE HOLDERS


SERRATED

21-8019	5"
21-8020	6"
21-8021	7"
21-8022	8"
21-8023	9"
21-8024	10¼"



TC OLSEN-HEGAR NEEDLE HOLDERS

21-8026	4¾"	Smooth	Delicate
21-8025	4¾"	Serrated	Delicate
21-8027	5½"	Serrated	
21-8028	6½"	Serrated	
21-8029	7½"	Serrated	




TC JACOBSON MICROVASCULAR NEEDLE HOLDERS

SMOOTH WITH LOCK

21-8053	7"	Straight
21-8253	7"	Curved
21-8055	8¼"	Straight
21-8255	8½"	Curved
21-8057	9"	Straight
21-8257	9"	Curved

WITHOUT LOCK

21-8054	7"	Straight
21-8254	7"	Curved
21-8058	8¼"	Straight
21-8258	8½"	Curved
21-8059	9"	Straight
21-8259	9"	Curved



TC JACOBSON MICROVASCULAR NEEDLE HOLDERS

SERRATED WITH LOCK

21-8250	7"	Straight
21-8251	7"	Curved
21-8052	8½"	Straight
21-8056	8¼"	Curved
21-8157	9"	Straight
21-8158	9"	Curved

WITHOUT LOCK

21-8150	7"	Straight
21-8151	7"	Curved
21-8152	8½"	Straight
21-8156	8¼"	Curved
21-8159	9"	Straight
21-8160	9"	Curved



TC CASTROVIEJO NEEDLE HOLDERS

STANDARD-WITH LOCK


- 21-8040 5½" Smooth, Straight
- 21-8046 5½" Serrated, Straight
- 21-8047 5½" Serrated, Curved



TC CASTROVIEJO NEEDLE HOLDERS

EXTRA DELICATE-WITHOUT LOCK


- 21-8044 5½" Smooth, Straight
- 21-8045 5½" Smooth, Curved



TC CASTROVIEJO NEEDLE HOLDERS

EXTRA DELICATE-WITH LOCK

- 21-8042 5½" Smooth, Straight
- 21-8043 5½" Smooth, Curved



TC MATHIEU NEEDLE HOLDERS

- 21-8070 5½"
- 21-8071 6¾" Narrow Jaw
- 21-8072 7½"

SKLARLITE THE LIGHT TOUCH

SKLARLITE XD EXTRA DELICATE

SKLARLITE™ AND SKLARLITE XD™ PROVIDE LIGHTNESS AND DURABILITY

In response to industry demands for lighter pattern instruments, Sklar has established the SklarLite line of premium operating room instruments. Through a unique manufacturing technique, we are able to provide a noticeably lighter instrument that maintains the strength and durability of a standard pattern. These extremely sensitive instruments give the surgeon maximum control and precision while comfortably withstanding the daily rigors of frequent use.


Sklarlite™ and Sklarlite XD™ premium instruments are our lightest and most sensitive instruments. These surprisingly strong instruments carry the same day-in and day-out reliability that you have come to expect from all Sklar instruments and are fully guaranteed for a lifetime of dependable use.

SKLARLITE THE LIGHT TOUCH

SKLARLITE™ DERF NEEDLE HOLDER 23-2124 4¾" Serrated

SKLARLITE XD EXTRA DELICATE

SKLARLITE XD™ DERF NEEDLE HOLDER 23-2224 4¾" Serrated



SKLARLITE THE LIGHT TOUCH


SKLARLITE™ WEBSTER NEEDLE HOLDERS

23-2125 5¼" Smooth
23-2120 5¼" Serrated

SKLARLITE XD EXTRA DELICATE

SKLARLITE XD™ WEBSTER NEEDLE HOLDERS

23-2225 5¼" Smooth
23-2220 5¼" Serrated



SKLARLITE THE LIGHT TOUCH


SKLARLITE™ HALSEY NEEDLE HOLDERS

23-2135 5¼" Smooth
23-2130 5¼" Serrated

SKLARLITE XD EXTRA DELICATE


SKLARLITE XD™ HALSEY NEEDLE HOLDERS

23-2235 5¼" Smooth
23-2230 5¼" Serrated



NEEDLE HOLDERS

Skarlite™



SKLARLITE
THE LIGHT TOUCH

**SKLARLITE™
BAUMGARTNER
NEEDLE HOLDER**

23-2132 5¼" Serrated

SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE XD™
BAUMGARTNER
NEEDLE HOLDER**

23-2232 5¼" Serrated



SKLARLITE
THE LIGHT TOUCH


**SKLARLITE™
CRILE-WOOD
NEEDLE HOLDER**

23-2134 6" Serrated

SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE XD™
CRILE-WOOD
NEEDLE HOLDER**

23-2234 6" Serrated



SKLARLITE
THE LIGHT TOUCH


**SKLARLITE™
MAYO-HEGAR
NEEDLE HOLDER**

23-2128 7" Serrated

SKLARLITE XD
EXTRA * DELICATE

**SKLARLITE XD™
MAYO-HEGAR
NEEDLE HOLDER**

23-2228 7" Serrated



SKLARLITE
THE LIGHT TOUCH


**SKLARLITE™
OLSEN-HEGAR
NEEDLE HOLDER**

23-2136 5½" Serrated

SKLARLITE XD
EXTRA * DELICATE


**SKLARLITE XD™
OLSEN-HEGAR
NEEDLE HOLDER**

23-2236 5½" Serrated



DERF NEEDLE HOLDER


20-1547 4¾" Serrated



WEBSTER NEEDLE HOLDERS

DELICATE


20-1350 5¼" Smooth
20-1352 5¼" Serrated



BROWN NEEDLE HOLDERS


ROUNDED JAWS

20-4852 5¼"
20-4867 6¾"




HALSEY NEEDLE HOLDERS

20-1750 5¼" Smooth
20-1751 5¼" Serrated




NEIVERT NEEDLE HOLDER

20-1760 5" One offset Ring



BAUMGARTNER NEEDLE HOLDER


20-1450 5¼" Serrated



CRILE-WOOD NEEDLE HOLDERS

SERRATED


20-1265	6"	Narrow Jaws
20-1270	7"	DeBakey Type
20-1280	8"	DeBakey Type
20-1290	9"	DeBakey Type
20-1295	10½"	DeBakey Type



MAYO-HEGAR NEEDLE HOLDERS

SERRATED


20-2052	5¼"	BROAD JAW
20-2060	6"	20-2270 7"
20-2070	7"	MASSON
20-2480	8"	20-2910 10½"
20-2490	9½"	20-2911 11½"



**OLSEN-HEGAR
NEEDLE HOLDERS**

SCISSORS COMBINATION


20-2748	4¾"	Smooth, Delicate
20-2747	4¾"	Serrated, Delicate
20-2755	5½"	Serrated
20-2765	6½"	Serrated



COLLIER NEEDLE HOLDER

FENESTRATED JAW

65-6550	5¼"	Serrated
---------	-----	----------



**JACOBSON MICRO NEEDLE
HOLDERS**

WITHOUT CATCH

98-3070	7"	Straight, Smooth
98-3085	8½"	Straight, Smooth

WITH CATCH

98-3170	7"	Straight, Smooth
98-3171	7"	Curved, Smooth



**CASTROVIEJO NEEDLE
HOLDERS**

WITH CATCH

65-6357	5½"	Straight
65-6358	5½"	Curved



**CASTROVIEJO
NEEDLE HOLDERS**


**STANDARD HANDLE, DELICATE, SMOOTH JAWS TO
HOLD 7-0 SUTURES**

WITHOUT CATCH

65-6350	5½"	Straight Jaw
65-6351	5½"	Curved Jaw


WITH CATCH

65-6352	5½"	Straight Jaw
65-6354	5½"	Curved Jaw



MATHIEU NEEDLE HOLDERS


20-3055	5½"	Serrated
20-3070	6¼"	Serrated
20-3075	7½"	Serrated



WEITLANER SELF-RETAINING RETRACTORS

BLUNT

22-8345	4½"	2x3 Teeth
22-8355	5½"	3x4 Teeth
42-1865	6½"	3x4 Teeth
42-1875	8"	3x4 Teeth
42-1895	9½"	3x4 Teeth




WEITLANER SELF-RETAINING RETRACTORS

SHARP

22-8245	4½"	2x3 Teeth
22-8255	5½"	3x4 Teeth
42-1965	6½"	3x4 Teeth
42-1975	8"	3x4 Teeth
42-1995	9½"	3x4 Teeth


SOLID BLADES

22-8760	5½"
---------	-----




GELPI RETRACTORS

22-8850	3½"
22-8853	5½"
22-8855	6¼"
60-6570	7"
60-6580	7½" With Stop




SCHINK METATARSAL SPREADER

97-5045	4½"
---------	-----




INGE RETRACTORS

42-1420	6½"
42-1425	10½"



ALM RETRACTORS

22-9027	2¾"	4x4 Teeth	Blunt
60-6340	4"	4x4 Teeth	Blunt
22-9029	2¾"	4x4 Teeth	Sharp
60-6342	4"	4x4 Teeth	Sharp




HOLZHEIMER SELF-RETAINING RETRACTORS

SHARP

22-8920 4" 4x4 Teeth

BLUNT

22-8940 4" 4x4 Teeth



VOLKMAN RETRACTORS

BLUNT, 4½"

60-4252 2 Prong
60-4253 3 Prong
60-4254 4 Prong
60-4256 6 Prong

SHARP, 4½"


60-4242 2 Prong
60-4243 3 Prong
60-4244 4 Prong
60-4246 6 Prong

BLUNT, 8½"

60-5502 2 Prong
60-5503 3 Prong
60-5504 4 Prong
60-5506 6 Prong

SHARP, 8½"

60-5602 2 Prong
60-5603 3 Prong
60-5604 4 Prong
60-5606 6 Prong



SENN RETRACTORS

DOUBLE END

60-6162 6¼" Sharp
60-6164 6¼" Blunt



RAGNELL RETRACTOR

DOUBLE END

22-8560 6"



DESMARRES LID RETRACTORS


65-1841 #1 12mm
65-1842 #2 14mm
65-1843 #3 16mm
65-1844 #4 18mm



ZALKIND (RIBBON) RETRACTORS

MALLEABLE

60-3594	1/4"x7 3/4"	60-3607	3/4"x13"
60-3596	1/2"x7 3/4"	60-3610	1"x13"
22-8067	5/8"x7 3/4"	60-3612	1 1/4"x13"
60-3598	3/4"x7 3/4"	60-3615	1 1/2"x13"
22-8075	1"x7 3/4"	60-3617	1 3/4"x13"
		60-3620	2"x13"
		60-3625	2 1/2"x13"



CRILE RETRACTOR

DOUBLE END

57-1445 4 1/2" Small



LUER RETRACTOR

DOUBLE END

22-8565 5 1/2"




MINI HOHMANN BONE ELEVATORS

40-7450	6 1/2"x6mm
40-7451	6 1/2"x8mm




HOHMANN BONE ELEVATORS

40-7447	8 1/2"	8mm
40-7448	9"	10mm
40-7452	9 1/2"	18mm
40-7453	9 1/2"	43mm
40-7454	10"	24mm
40-7455	10"	20mm
40-7456	10"	22mm
40-7458	10 1/2"	65mm



TRACHEAL RETRACTORS


6 1/2"		
78-1001	Sharp	1 Prong
78-1002	Sharp	2 Prongs
78-1003	Sharp	3 Prongs
78-1101	Blunt	1 Prong
78-1102	Blunt	2 Prongs
78-1103	Blunt	3 Prongs



PLASTIC SURGERY RETRACTORS


6" FLEXIBLE SHANK

22-9144	Sharp	1 Prong
22-9132	Sharp	2 Prongs
22-9133	Sharp	3 Prongs
22-9141	Blunt	1 Prong
22-9142	Blunt	2 Prongs
22-9143	Blunt	3 Prongs



JOSEPH SKIN HOOKS SINGLE END, SHARP PRONGS, 6¼"


47-2500	1 Prong
47-2502	2 Prongs 2mm
47-2505	2 Prongs 5mm
47-2507	2 Prongs 7mm
47-2510	2 Prongs 10mm
47-2512	2 Prongs 12mm



GUTHRIE DOUBLE HOOKS


SHARP PRONGS

47-2651	2 Small
47-2650	2 Large



FREER SKIN HOOKS DOUBLE PRONG, 2.5MM


47-2760	6"	Sharp
47-2758	6"	Blunt



BLAIR (ROLLET) RETRACTORS


4 PRONG-DELICATE

47-2655	5½"	Sharp
47-2752	5½"	Blunt



GILLIES SKIN HOOKS


47-2570	7"	3mm
47-2572	7"	4mm



LAHEY RETRACTORS

4¾" SHARP

78-1201	1 Prong
78-1202	2 Prongs
78-1203	3 Prongs



TYRELL IRIS HOOKS

65-3205	5"	Blunt
65-3206	5"	Sharp



FRAZIER TISSUE HOOKS

EXTRA DELICATE

47-2555 5" Sharp
47-2556 5" Blunt



SKIN HOOK


SHARP

47-2550 4¾" 2mm



BERGMAN HOOK

47-2302 5½" Sharp




KILNER ALAR HOOKS

47-2520 3¼" Sharp
47-2575 5¾" Sharp



O'CONNOR HOOK

65-3265 5" Sharp



RAMELET PHLEBECTOMY HOOKS

- 06-4170 Right, Fine
- 06-4173 Left, Fine
- 06-4175 Set of Two, 1 Fine, 1 Medium, Left
- 06-4171 Right, Medium
- 06-4174 Left, Medium
- 06-4172 Set of Two, 1 Fine, 1 Medium, Right




WIENER CORNEAL HOOKS

- 65-3250 5" Sharp
- 65-3248 5" Blunt




MCGANNON IRIS RETRACTOR

- 65-3258 5 1/4"




JAEGER LID PLATE

- 65-1835 4 1/4" Stainless




TC SKLAR STERNAL WIRE TWISTERS

21-9060 5"
21-9065 7"




SKLAR EDGE™ TC WIRE CUTTING SCISSOR

16-2800 4¾" Serrated Angled




TC PIN CUTTERS SIDE CUT DOUBLE ACTION

40-1023 7" Straight
40-1024 8½" Angular




TC PIN CUTTER DOUBLE ACTION

Cuts to .062" Wire 1.6mm
40-1045 7" End Cut-Side Cut




TC WIRE AND PIN CUTTER

40-1018 8"



TC WIRE PULLER

40-1730 8"




TC WIRE PULLER

DOUBLE ACTION


40-1735 7" 4mm (A)

40-1738 7" 6mm (B)




TC PIN PULLER

40-1740 5½"



**TC RING HANDLE WIRE
EXTRACTOR**

40-1720 7"



**SUTURE WIRE CUTTING
SCISSORS**

ANGULAR SERRATED


24-2342 4¾"



**DOUBLE ACTION PIN
CUTTER**

STAINLESS


40-1715 15" Cuts up to ⅝"




DIAMOND PIN CUTTER

CHROME


40-1762 6¼"




SKLAR PIN CUTTER
END CUTTING
40-1048 7"




IMPROVED SKLAR FORCEPS
FOR CUTTING AND SHAPING WIRE
75-9550




FLAT NOSE PLIER
40-1778 5½"




NEEDLE NOSE PLIER
40-1775 5½"



DIAMOND PIN CUTTER
40-1762 6¼" Cuts up to 5/64"




SNAP LOCK PIN/WIRE EXTRACTOR
40-1725



SKLAR NAIL DRILL

THREE DRILLS IN HANDLE


- 06-4160 Complete
- 06-4161 Set of Drill Points Only



SKLAR BONE DRILL-OFFSET HANDLE

WITH JACOBS CHUCK

- 40-1465 For Twist Drills and Kirschner Wires
- 40-1478 Chuck Key Only



BUNNELL DRILL-CANNULATED WITH JACOBS CHUCK

- 40-1475 For Twist Drills and Kirschner Wires
- 40-1478 Chuck Key Only



KIRSCHNER WIRES

TWO TROCAR ENDS


- | | | | |
|---------|------------|---------|------------|
| 40-1565 | 4" x .028 | 40-1580 | 9" x .035 |
| 40-1566 | 4" x .035 | 40-1582 | 9" x .045 |
| 40-1567 | 4" x .045 | 40-1584 | 9" x .062 |
| 40-1568 | 4" x .054 | 40-1535 | 12" x .035 |
| 40-1569 | 4" x .062 | 40-1545 | 12" x .045 |
| 40-1571 | 7½" x .045 | 40-1562 | 12" x .062 |
| 40-1572 | 7½" x .062 | | |



STEINMANN PINS

TROCAR POINT

- | | | | |
|---------|------------|---------|------------|
| 40-1590 | 9" x 2mm | 40-1594 | 9" x 3.5mm |
| 40-1591 | 9" x 2.4mm | 40-1595 | 9" x 4.0mm |
| 40-1592 | 9" x 2.8mm | 40-1596 | 9" x 4.5mm |
| 40-1593 | 9" x 3.2mm | | |



OSCILLATING PLASTER SAWS

11-2100	110V Complete	
11-2101	220V Complete	
11-2105	Blade Only	50mm
11-2110	Blade Only	65mm
11-2112	Blade Only	65mm, For use on Synthetic Plaster
11-2115	Wrenches only	Set of 2



WALTON PLASTER SPREADER

3 BLADES
11-1990 9"




HENNIG PLASTER SPREADER

2 BLADES
11-1910 11"



BEESON PLASTER SPREADER

2 BLADES
11-1812 12"




**WOLFE-BOEHLER PLASTER
CAST BREAKER**

11-2095 10"



PEDIATRIC CAST BREAKER

11-2098 7" (18 cm) Mason




ENGEL PLASTER SAW

11-1525 Chrome Plated



WIGMORE PLASTER SAW

11-1530 7½"




REINER PLASTER KNIFE

METAL HANDLE
11-1550 7" 1½" Blade



MURPHY PLASTER KNIFE


WOOD HANDLE
11-1540



SEUTIN PLASTER SHEARS

HEAVY PLIER HANDLE


11-1490 9" Serrated



BRUNS PLASTER SHEARS


SERRATED

11-1295 9¼"
16-2820 9½" TC



HERCULES UTILITY SCISSOR


11-1297 7½" Serrated



COMMANDER UTILITY SHEARS

HEAVY DUTY WITH WIRE CUTTING NOTCH


11-1298 8½"



BRUNS PLASTER SHEARS

HEAVY PATTERN


11-1380 8"
11-1390 9"



STILLE PLASTER SHEARS

11-1610 10"


11-1615 15"



**TAYLOR PERCUSSION
HAMMERS**

SOLID HANDLE

06-3170 Small 7"
06-3175 Large 7½"



**RABINER NEUROLOGICAL
HAMMER**

WITH PIN AND BRUSH


06-3290 Chrome



**BUCK NEUROLOGICAL
HAMMER**


WITH PIN AND BRUSH

06-3370 6¾"



ALUMINUM ALLOY TUNING FORKS

67-7064	C-64	67-7524	C-1024
67-7128	C-128	67-7548	C-2048
67-7256	C-256	67-7596	C-4096
67-7512	C-512		




TUNING FORK SET

ALUMINUM ALLOY

Set of 5
67-7605 Including C-128 C-1024
C-256 C-2048
C-512



**WARTENBERG NEUROLOGICAL
PINWHEEL**
CHROME PLATED
06-3375



TAPE MEASURE RETRACTABLE
06-2852 60"



STEEL RULER
06-2860 6"
06-2862 8"
06-2861 12"
06-2863 20"


Cleanwheel™


STERILE • DISPOSABLE
Neurological Pinwheel for
Sensory Evaluation

Product Numbers

96-7905 Box of 36

96-7907 Case of 4 Boxes, 144 Ct.





TROUTMAN BLADE BREAKER

65-2960 4"



CASTROVIEJO BLADE BREAKER


65-2950 4"



SWISS MODEL BLADE BREAKER

4 3/4"

65-2958 15x2.5mm Jaws




KEYES FINGER RING CUTTERS

STAINLESS STEEL

06-4150 Ring Cutter Complete

06-4151 Blades w/Brass Bushings Attached for above Ring Cutter



SKLAR FINGER RING CUTTERS

STAINLESS STEEL


06-4152 6"

06-4153 Blade Only

SCALPEL HANDLES

PRODUCT NUMBER FEATURES

06-2910 #3 with metric scale



PRODUCT NUMBER FEATURES

06-2903 #3 fits Blade #10, 11, 12, 15

06-2915 #3L

06-2916 #3XL



PRODUCT NUMBER FEATURES

06-2904 #4 fits Blade #20, 21, 22, 23

06-2920 #4L



PRODUCT NUMBER FEATURES

06-2907 #7 fits Blade #10, 11, 12, 15



PRODUCT NUMBER FEATURES

06-2911 #9 fits Blade #10, 11, 12, 15



PRODUCT NUMBER FEATURES

06-2912 #8 fits Blade #60



PRODUCT NUMBER FEATURES

06-2925 Large Lab Type Handle #3

06-2930 Large Lab Type Handle #4



PRODUCT NUMBER FEATURES


06-2908 Round Scalpel Handle, Straight
#3R, fits Blade #10, 11, 12, 15

06-2909 Round Scalpel Handle, Angled
Fits Blade #10, 11, 12, 15



PRODUCT NUMBER FEATURES

06-3090 Blade Remover





BLADES, STAINLESS STEEL

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3100	#10	100/box
06-3101	#11	100/box
06-3102	#12	100/box
06-3103	#15	100/box
06-3108	#15C	100/box
06-3104	#20	100/box
06-3105	#21	100/box
06-3106	#22	100/box
06-3107	#23	100/box

BLADES, CARBON STEEL

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3010	#10	100/box
06-3011	#11	100/box
06-3012	#12	100/box
06-3013	#12D	100/box
06-3015	#15	100/box
06-3016	#15C	100/box
06-3020	#20	100/box
06-3021	#21	100/box
06-3022	#22	100/box
06-3023	#23	100/box

SCALPELS, DISPOSABLE

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3110	#10	10/box
06-3111	#11	10/box
06-3112	#12	10/box
06-3115	#15	10/box
06-3120	#20	10/box
06-3121	#21	10/box
06-3122	#22	10/box
06-3123	#23	10/box



SKLARSAFE™ SAFETY SCALPELS



LATEX-FREE

- Protective, transparent sheath can be opened, closed, and locked with a single hand.
- Precision stainless steel blade made in England.
- Sturdy, “Sure-Grip” handle.
- Audible “click” when sheath advanced to locked position.
- OSHA Compliant.
- Individually sealed in sturdy film and Tyvek® pack.

STERILE

PRODUCT NUMBER FEATURES QUANTITIES

06-3094 #10 10/box

06-3095 #11 10/box

06-3096 #15 10/box

NON-STERILE

06-3130 #10 500/cs

06-3132 #11 500/cs

06-3134 #15 500/cs

MINI BLADES AND SCALPELS

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
97-1822	#61	25/box
97-1821	#62	25/box
97-1820	#63	25/box
97-1826	#64	25/box
97-1825	#65	25/box
97-1824	#67	25/box
97-1823	#69	25/box



MINI PODIATRY HANDLE

NON-STERILE


PRODUCT NUMBER	FEATURES	QUANTITIES
97-1450	#3K	1 Each



SPECIALTY BLADES

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
97-1400	#318, Curved Edge	12/box
97-1405	#312, Curved Edge	12/box
97-1410	#313, Curved Edge	12/box
97-1420	#316, Straight Edge	12/box



SELF-LOCKING HANDLE

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
97-1445	#6	1 Each



CARBON P-M BLADE

NON-STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3024	#60	50/box



STITCH CUTTER

STERILE

PRODUCT NUMBER	FEATURES	QUANTITIES
06-3030	Stainless Steel	100/box





SKLAR TRU-PUNCH™ DISPOSABLE BIOPSY PUNCH

The SKLAR Tru-punch has a seamless steel tip carefully honed for a razor sharp incision and perfect sample. The ribbed and hollow handle permits tissue to be easily dislodged.

STERILE DISPOSABLE BIOPSY PUNCH

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-0192	1mm	25/bx
96-1097	1.5mm	25/bx
96-1140	2mm	25/bx
96-0102	2.5mm	25/bx
96-1142	3mm	25/bx
96-1144	3.5mm	25/bx
96-1146	4mm	25/bx
96-1148	5mm	25/bx
96-1150	6mm	25/bx
96-1152	8mm	25/bx
96-1160	Assorted - 5 of each: 2, 3, 3.5, 4 & 5mm	25/bx

PRODUCT NUMBER	DESCRIPTION	QUANTITIES
96-1092	1mm	50/bx
96-1098	1.5mm	50/bx
96-1100	2mm	50/bx
96-1102	2.5mm	50/bx
96-1105	3mm	50/bx
96-1110	3.5mm	50/bx
96-1115	4mm	50/bx
96-1120	5mm	50/bx
96-1125	6mm	50/bx
96-1130	8mm	50/bx
96-0160	Assorted - 5 of each: 2, 3, 3.5, 4 & 5mm	50/bx


SKLAR DISPOSABLE DERMAL CURETTES SURGICAL STEEL TIPS



STERILE PRODUCT NUMBER	FEATURES	QUANTITIES
96-1190	2 mm	25/bx
96-1192	3 mm	25/bx
96-1194	4 mm	25/bx
96-1196	5 mm	25/bx
96-1198	7 mm	25/bx
96-1199	Assorted, 5 each	25/bx

ECONO™ STERILE SINGLE-USE INSTRUMENTS

- Greatly Reduces HAIs & Cross Contamination
- Eliminates Costly Processing, Cleaning & Sterilization



ECONO™ INSTRUMENTS
A STERILE, SINGLE-USE INSTRUMENT
HALSTEAD MOSQUITO FORCEPS, STRAIGHT, 5"
LOT # 116999
* INSTRUMENT IS STERILE UNLESS PACKAGE HAS BEEN OPENED OR DAMAGED
* FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN
ECONO™ BRAND IMPORTED FROM PARSIPpany, NJ
SKLAR INSTRUMENTS • 889 SOUTH MATLACK ST. • WEST CHESTER, PA 19380
PH: (610) 430-3200 • FAX: (610) 429-0500
E-Mail: surg@sklarcorp.com • Website: <http://www.sklarcorp.com>

**CONTACT YOUR SKLAR ACCOUNT
MANAGER FOR DERMAL AND PODIATRY
STERILE DISPOSABLE INSTRUMENTS**

Available in boxes of 25, 50 and 100

98-1058 BASIC DERMAL INSTRUMENT SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
06-4120	Shamberg Comedone Extractor	Pinched End	1
06-4125	Unna Comedone Extractor	5¾"	1
11-1045	Lister Bandage Scissors	4½"	1
11-1055	Lister Bandage Scissors	5½"	1
15-1555	Mayo Dissecting Scissors	Straight, 5½"	1
15-2555	Mayo Dissecting Scissors	Curved, 5½"	1
17-1235	Hartmann Mosquito Forceps	Straight, Serrated, 3½"	1
17-1335	Hartmann Mosquito Forceps	Curved, Serrated, 3½"	1
17-1450	Halsted Mosquito Forceps	Straight, 5"	1
17-1550	Halsted Mosquito Forceps	Curved, 5"	1
19-1055	Dressing Forceps	Serrated, 5½"	1
19-1255	Tissue Forceps	1x2 Teeth, 5½"	1
19-3035	Fine Point Splinter Forceps	3½"	1
19-3045	Fine Point Splinter Forceps	4½"	1
20-1350	Webster Needle Holder	Smooth, 5¼"	1
20-1352	Webster Needle Holder	Serrated, 5¼"	1
20-1547	Derf Needle Holder	Serrated, 4¾"	1
20-1750	Halsey Needle Holder	Smooth, 5"	1
20-1751	Halsey Needle Holder	Serrated, 5"	1
20-2747	Olsen-Hegar Needle Holder	4¾"	1
20-2755	Olsen-Hegar Needle Holder	5½"	1
30-1945	Probe	With Eye, 4½"	1
30-1955	Probe	With Eye, 5½"	1
36-2060	Allis Tissue Forceps	3x4 Teeth, 6"	1
36-2275	Allis Tissue Forceps	5x6 Teeth, 7½"	1
47-1145	Iris Scissors	Straight, Sharp/Sharp, 4½"	1
47-1245	Iris Scissors	Curved, Sharp/Sharp, 4½"	1
47-1246	Iris Scissors	Straight, Sharp/Sharp, Delicate, 4½"	1
47-1248	Iris Scissors	Curved, Sharp/Sharp, Delicate, 4½"	1
47-1360	Metzenbaum Scissors	Curved, 5¾"	1
47-2047	Adson Tissue Forceps	1x2 Teeth, Delicate, 4¾"	1
47-2147	Adson Forceps	Serrated, Delicate, 4¾"	1
47-2935	Backhaus Towel Clamp	3½"	1
47-2955	Backhaus Towel Clamp	5¼"	1
50-3047	Adson Tissue Forceps	1x2 Teeth, 4¾"	1
50-3147	Adson Dressing Forceps	Serrated, 4¾"	1
64-1242	Stevens Tenotomy Scissors	Straight, Blunt/Blunt, 4½"	1

98-1058 BASIC DERMAL INSTRUMENT SET (CONTINUED)

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
64-1342	Stevens Tenotomy Scissors	Curved, Blunt/Blunt, 4½"	1
67-2502	Buck Ear Curette	#2, Angled, Sharp	1
67-2517	Buck Ear Curette	#2, Straight, Sharp	1
74-1370	Ballenger Sponge Forceps	Straight, Smooth, 7"	1
74-1670	Ballenger Sponge Forceps	Curved, Serrated, 7"	1
75-5460	Metzenbaum Scissors	Straight, 5¾"	1

98-1015 BASIC DERMAL EXCISION SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
06-2903	Scalpel Handle	#3	1
06-3011	Surgical Blades	#11, Box of 100, Carbon Steel	1
06-3015	Surgical Blades	#15, Box of 100, Carbon Steel	1
20-1750	Halsey Needle Holder	Smooth, 5"	1
22-2545	Metzenbaum Scissors	Curved, 4¾"	1
47-1145	Iris Scissors	Straight, Sharp/Sharp, 4½"	1
47-1245	Iris Scissors	Curved, Sharp/Sharp, 4½"	1
47-2047	Adson Tissue Forceps	1x2 Teeth, Delicate, 4¾"	1
47-2147	Adson Forceps	Serrated, Delicate, 4¾"	1

98-1645 MAJOR PODIATRY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
06-2903	Scalpel Handle	#3	1
11-1055	Lister Bandage Scissorss	5 1/2"	1
14-1055	Operating Scissorss	Straight, Sharp/Blunt, 5 1/2"	1
17-1450	Halsted Mosquito Forceps	Straight, 5"	6
17-1550	Halsted Mosquito Forceps	Curved, 5"	6
17-2155	Kelly Forceps	Curved, 5 1/2"	4
18-2155	Allis Tissue Forceps	Delicate, 4x5 Teeth, 5"	2
19-1255	Tissue Forceps	1x2 Teeth, 5 1/2"	2
20-1547	Derf Needle Holder	Serrated, 4 3/4"	1
20-2060	Mayo-Hegar Needle Holder	Serrated, 6"	1
22-8067	Zalkind Ribbon Retractor	5/8"x7 3/4"	2
22-8245	Weitlaner Retractor	Sharp, 2x3 Teeth, 4 1/2"	1
24-2340	Crown Scissorss	Straight, Serrated, 4"	1
40-3040	Liston Bone Forceps	Straight with Sharp Tips, 5 1/2"	1
40-4560	Boehler Bone Cutting Forceps	Straight, Double Action, 6"	1
40-5875	Lucae Bone Mallet	Small, 7 1/2"x1"	1
40-6552	Lambotte Ostetomes	Straight, 5", Set of 9	1
40-7250	Volkman Bone Curette	Double End, 5"	1
40-9688	Fomon Rasp	4 Sided, Convex/Flat, 8 1/2"	1
40-9825	Joseph Rasp	Straight, Fine Cross Serrations, 6 1/4"	1
47-1245	Iris Scissorss	Curved, Sharp/Sharp, 4 1/2"	1
47-1447	Adson Brown Forceps	9x9 Teeth, Side Grasping, 4 3/4"	1
47-2935	Backhaus Towel Clamp	3 1/2"	3
47-2955	Backhaus Towel Clamp	5 1/4"	1
50-3047	Adson Tissue Forceps	1x2 Teeth, 4 3/4"	3
60-6162	Senn Retractors	Double-End, Sharp, 3 Prongs, 6 1/4"	2
65-3206	Tyrell Iris Hook	Sharp, 5"	3
87-2195	Foerster Sponge Forceps	Straight, Serrated, 9 1/2"	3
97-0623	Locke Elevator	Narrow, 3x15mm	1
97-1318	Platypus Nail Forceps	6mm, 5 1/2"	1

98-1010 BASIC PODIATRY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
06-2903	Scalpel Handle	#3	1
06-3011	Surgical Blades	#11, Box of 100	1
06-3015	Surgical Blades	#15, Box of 100	1
11-1055	Lister Bandage Scissorss	5 1/2"	1
11-1255	Knowles Bandage Scissorss	5 1/2"	1
13-1055	Operating Scissorss	Straight, Sharp/Sharp, 5 1/2"	1
15-1355	Stitch Scissorss	5 1/2"	1
16-2800	TC Wire Cutting Scissorss	Angled, 4 3/4"	1
17-1450	Halsted Mosquito Forceps	Straight, 5"	3
17-1550	Halsted Mosquito Forceps	Curved, 5"	3
19-1055	Dressing Forceps	Serrated, 5 1/2"	2
19-1255	Tissue Forceps	1x2 Teeth, 5 1/2"	1
19-3045	Fine Point Splinter Forceps	4 1/2"	1
21-8002	TC Webster Needle Holder	Smooth, 5 1/4"	1
21-8027	TC Olsen Hegar Needle Holder	Serrated, 5 1/2"	1
23-1102	Sklarlite™ Operating Scissorss	Straight, Sharp/Sharp, 5 1/2"	1
23-1106	Sklarlite™ Operating Scissorss	Straight, Sharp/Blunt, 5 1/2"	1
23-1189	Sklarlite™ Spencer Stitch Scissorss	4 1/2"	1
23-2110	Sklarlite™ Crile Forceps	Straight, 5 1/2"	3
23-2113	Sklarlite™ Crile Forceps	Curved, 5 1/2"	3
23-2602	Sklarlite™ Adson Tissue Forceps	1x2 Teeth, 4 3/4"	2
40-9825	Joseph Rasp	Straight, Fine Cross Serrations, 6 1/2"	1
47-2380	Freer Elevator	Double-End, Sharp/Blunt, 7 1/4"	1
47-2935	Backhaus Towel Clamp	3 1/2"	6
50-1355	Dandy Scalp Hemostatic Forceps	5 1/2"	1
50-2006	Frazier Suction Tube	6 French	2
50-3147	Adson Dressing Forceps	Serrated, 4 3/4"	2
97-0335	Tissue Nipper	Convex Jaw, 9mm, 4"	1
97-0532	Curette/Excavator	#58 1-2, Double End	1
97-0536	Curette/Excavator	#58 1-3, Double End	1
97-0542	Spatula and Packer	#91, Double End	1
97-0545	Corn Chisel	Curved Edge, Spoon	1
97-0554	Ingrown Nail Shaver	Fenestated	1
97-1150	Ingrown Nail Splitter	English Anvil, 5"	1
97-1318	Platypus Nail Forceps	6mm, 5 1/2"	1

98-1650 GENERAL PODIATRY SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
06-2903	Scalpel Handle	#3	1
06-2904	Scalpel Handle	#4	1
11-1055	Lister Bandage Scissorss	5 1/2"	1
17-1450	Halsted Mosquito Forceps	Straight, 5"	1
17-1550	Halsted Mosquito Forceps	Curved, 5"	1
23-9259	Sklarlite Littauer Cutting Forceps	Straight Back, 5"	1
47-1245	Iris Scissorss	Curved, Sharp/Sharp, 4 1/2"	1
47-2955	Backhaus Towel Clamps	5 1/4"	2
50-3147	Adson Dressing Forceps	Serrated, 4 3/4"	1
97-0474	Nail Nipper	Concave, Double Spring with Lock, 5 1/2"	1
97-0536	Curette/Excavator	#58 1-3, Double End	1
97-0542	Spatula and Packer	#91, Double End	1
97-0545	Corn Chisel	Curved Edge, Spoon	1
97-1056	Nail Nipper	Double Spring with Lock, 5 1/2"	1
97-1150	Ingrown Nail Splitter	English Anvil, 5"	1
97-1262	Littauer Cutting Forceps	6"	1
97-1351	Moleskin Scissorss	Curved, Serrated, 7 1/2"	1

98-1655 SMALL BONE SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
06-2860	Ruler	Steel, 6"	1
06-2903	Scalpel Handle	#3	2
10-1675	Medicine Cup	2oz., Graduated Measure	1
15-1567	Mayo Dissecting Scissorss	Straight, 6 3/4"	1
17-1550	Halsted Mosquito Forceps	Curved, 5"	6
17-2862	Rochester-Ochsner Forceps	Straight, 1x2 Teeth, 6 1/4"	2
17-3155	Crile Forceps	Curved, 5 1/2"	2
18-2162	Allis Tissue Forceps	4x5 Teeth, 6"	2
19-1255	Tissue Forceps	1x2 Teeth, 5 1/2"	2
20-1265	Baby Crile-Wood Needle Holder	Serrated, 6"	2
20-1270	Crile-Wood Needle Holder	7"	1
22-8245	Weitlaner Retractor	Sharp, 2x3 Teeth, 4 1/2"	1
40-2580	Kern Bone Forceps	with Ratchet, 5 3/4"	2
40-2585	Kern Bone Forceps	with Ratchet, 8 1/2"	2
40-2670	Lewin Bone Forceps	Serrated, 7"	2
40-2921	Lowman Clamp	1x2 Prongs, 4 3/4"x 3/4" Jaw	2

98-1655 SMALL BONE SET (CONTINUED)

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
40-2922	Lowman Clamp	1x2 Prongs, 7 1/4" x 1 1/8" Jaw	2
40-3055	Liston Bone Forceps	Straight, 5 1/2"	1
40-4560	Boehler Bone Cutting Forceps	Straight, Double Action, 6"	1
40-5765	Gerzog Bone Mallet	7 1/2"	1
40-6953	Hoke Osteotome	Straight, 3mm, 5 1/2"	1
40-6954	Hoke Osteotome	Straight, 4mm, 5 1/2"	1
40-6956	Hoke Osteotome	Straight, 6mm, 5 1/2"	1
40-6958	Hoke Osteotome	Straight, 8mm, 5 1/2"	1
40-6970	Hoke Osteotome	Straight, 10mm, 5 1/2"	1
40-6973	Hoke Osteotome	Straight, 13mm, 5 1/2"	1
40-7400	Key Elevator	1/4" x 7 1/2"	1
40-7402	Key Elevator	1/2" x 7 1/2"	1
40-7450	Mini Hohmann Bone Elevator	6mm, 6 1/2"	1
40-7451	Mini Hohmann Bone Elevator	8mm, 6 1/2"	1
40-7452	Mini Hohmann Bone Elevator	18mm, 9 1/2"	1
47-1047	Wagner Scissorss	Straight, Sharp/Sharp, 4 3/4"	1
47-1058	Wagner Scissorss	Curved, Sharp/Sharp, 4 1/2"	1
47-1360	Metzenbaum Scissorss	Curved, 5 3/4"	1
47-1447	Adson Brown Forceps	9x9 Teeth, Side Grasping, 4 3/4"	2
47-2380	Freer Elevator	Double-End, Sharp/Blunt, 7 1/4"	1
47-2502	Joseph Hook	2 Sharp Prongs, 2mm, 6 1/4"	2
47-2935	Backhaus Towel Clamp	3 1/2"	2
49-1355	Wax Spatula	Double End, #7	1
50-2009	Frazier Suction Tube	9 French	1
50-2011	Frazier Suction Tube	11 French	1
50-3047	Adson Tissue Forceps	1x2 Teeth, 4 3/4"	2
60-1085	US Army Retractors	Set of 2, Double End, 8 1/4"	2
60-5822	Volkman Retractor	Two Prongs, Sharp, 8 1/2"	2
60-6162	Senn Retractors	Double-End, Sharp, 3 Prongs, 6 1/4"	2

98-1660 TOENAIL REMOVAL SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
06-2903	Scalpel Handle	#3	1
17-1450	Halsted Mosquito Forceps	Straight, 5"	1
47-1447	Adson Brown Forceps	9x9 Teeth, Side Grasping, 4 3/4"	1
47-2385	McKenty Septum Elevator	Curved, 4mm, 6"	1
97-0532	Curette/Excavator	#58 1-2, Double End	1
97-0542	Spatula and Packer	#91, Double End	1
97-1150	Ingrown Nail Splitter	English Anvil, 5"	1
97-1242	Tissue Nipper	Convex, 4 1/2"	1
97-1315	Platypus Nail Forceps	4mm, 5 1/2"	1
97-1318	Platypus Nail Forceps	6mm, 5 1/2"	1
97-1450	Scalpel Handle	#3K For Miniature Edge Blades	1

98-1665 TOE PACK SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
14-1055	Operating Scissorss	Straight, Sharp/Blunt, 5 1/2"	1
15-1555	Mayo Dissecting Scissorss	Straight 5 1/2"	1
17-1450	Halsted Mosquito Forceps	Straight, 5"	1
17-1550	Halsted Mosquito Forceps	Curved, 5"	2
19-1255	Tissue Forceps	1x2 Teeth, 5 1/2"	1
19-1260	Tissue Forceps	1x2 Teeth, 6"	2
20-1547	Derf Needle Holder	Serrated, 4 3/4"	2
22-2652	Ragnell Dissecting Scissorss	Curved, 5"	1
40-4560	Boehler Bone Cutting Forceps	Straight, Double Action, 6"	1
40-9825	Joseph Rasp	Straight, Fine Cross Serrations, 6 1/4"	2
47-1360	Metzenbaum-Lahey Scissorss	Curved, 5 3/4"	1
47-1447	Adson Brown Forceps	9x9 Teeth, Side Grasping, 4 3/4"	2
47-2380	Freer Elevator	Double-End, Sharp/Blunt, 7 1/4"	1
47-2955	Backhaus Towel Clamp	5 1/4"	2
60-6162	Senn Retractor	Double-End, Sharp, 3 Prongs, 6 1/4"	2
65-3205	Tyrell Iris Hook	Blunt, 5"	2
67-5990	Blumenthal Rongeur	30 Degrees, Angled, 6"	1
87-2195	Foerster Sponge Forceps	Straight, Serrated, 9 1/2"	2
97-1450	Scalpel Handle	#3K For Miniature Edge Blades	1

98-1670 BUNION SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
06-2903	Scalpel Handle	#3	2
14-1055	Operating Scissorss	Straight, Sharp/Blunt, 5 1/2"	1
15-1555	Mayo Dissecting Scissorss	Straight 5 1/2"	1
17-1450	Halsted Mosquito Forceps	Straight, 5"	10
17-2855	Rochester-Ochsner (Kocher) Forceps	Straight, 1x2 Teeth, 5 1/2"	2
19-1255	Tissue Forceps	1x2 Teeth, 5 1/2"	2
19-1260	Tissue Forceps	1x2 Teeth, 6"	1
20-1547	Derf Needle Holder	Serrated, 4 3/4"	2
22-2760	Reynolds Scissorss	Curved, Narrow Tip, 6"	1
22-8067	Zalkind Ribbon Retractor	5/8"x7 3/4"	2
22-8245	Weitlaner Retractor	Sharp, 2x3 Teeth, 4 1/2"	1
40-4560	Boehler Bone Cutting Forceps	Straight, Double Action, 6"	1
40-9825	Joseph Rasp	Straight, Fine Cross Serrations, 6 1/2"	1
47-1360	Metzenbaum Scissorss	Curved, 5 3/4"	1
47-1447	Adson Brown Forceps	9x9 Teeth, Side Grasping, 4 3/4"	2
47-2380	Freer Elevator	Double-End, Sharp/Blunt, 7 1/4"	1
47-2955	Backhaus Towel Clamp	5 1/4"	2
60-6162	Senn Retractors	Double-End, Sharp, 3 Prongs, 6 1/4"	2
65-3205	Tyrell Iris Hook	Blunt, 5"	2
67-5990	Blumenthal Rongeur	30 Degrees, Angled, 6"	1
87-2195	Foerster Sponge Forceps	Straight, Serrated, 9 1/2"	2
97-1450	Scalpel Handle	#3K For Miniature Edge Blades	1

98-1675 SUTURE REMOVAL SET

PRODUCT NUMBER	DESCRIPTION	FEATURES	QUANTITIES
19-1055	Dressing Forceps	Serrated, 5 1/2"	1
19-1255	Tissue Forceps	1x2 Teeth, 5 1/2"	1
15-1355	Stitch Scissorss	5 1/2"	1

Sklar is pleased to offer a full line of Care & Cleaning agents to disinfect, decontaminate, sterilize, lubricate and protect your stainless steel surgical instruments and equipment. For expert care of your surgical instruments at every stage of the cleaning process, trust Sklar to deliver the right cleaning agent for you.



3 EASY STEPS FOR THE CARE AND CLEANING OF YOUR INSTRUMENTS

1 DECONTAMINATION 2 MINUTE KILL TIME

Sklar Soak™ and Sklar Disinfectant are intended to safely decontaminate soiled surgical instruments immediately after use. These ready to use, non-toxic formulas are safe for use on medical and surgical equipment including stainless steel, chrome, glass, plastic, Plexiglas, baked on painted and other non-porous surfaces. Effective against H1N1, HIV-1, HSV, VRE, Rotaviruses, Hep B & C as well as TB, MRSA and A-Baumannii.



2 CLEANING DETERGENTS AND ENZYMATICS

Choose from a wide selection of Sklar cleaning agents. First use Sklar Spray-Zyme™ foam pre-cleaner to coat instruments until ready to process. Next use Sklar Enzymatic™, a fast acting detergent specifically formulated for soaking/cleaning instruments heavily laden with organic debris. Finally, use Sklar Kleen™ as your neutral pH, non-corrosive detergent for manual, ultrasonic or automatic cleaning. Sklar detergents and enzymatics are effective in the cleaning process of surgical instruments, as well as metals, glassware, plastics and rubber.



3 LUBRICATION SKLAR LUBRICANT

To guarantee the extended life and performance of your surgical instruments, lubrication should be a standard part of your instrument processing regimen. Easy to use Sklar Instru-Guard™ products improve function and contain rust inhibitors that protect against corrosion, rust and stains. Following lubrication, instruments should be sterilized according to the instrument manufacturer's instructions and/or stored until use.



SKLAR-CIDE™

HIGH LEVEL DISINFECTING/STERILIZING SOLUTION

Sporicidal • Bactericidal • Virucidal • Fungicidal

SKLAR-CIDE provides a safe and efficient means for the high-level disinfection or sterilization of medical, ophthalmic, dental and surgical devices including endoscopes and lensed instruments. It may be used as a standard soak solution or in automatic reprocessing units set at room temperature. Unlike other solutions, there is no need to increase the solution temperature to provide high-level disinfection in 25 minutes.

SKLAR-CIDE is a surfactant-free, activated glutaraldehyde liquid chemical sterilant and high-level disinfectant. Its active ingredient is glutaraldehyde in a nominal concentration of 3%. When the powdered activator is added, the solution turns green in color to indicate its readiness for use. In addition to the active ingredient, this ready-to-use solution contains buffers, wetting agents, and a rust and corrosion inhibitor and it has a pleasant fragrance.

SKLAR-CIDE produces high-level disinfection at room temperature in as little time as 25 minutes, and complete sterilization in 10 hours at room temperature.

- Reusable for 28 days instead of 14, safer and saves time and money.
- Surfactant free - Will not damage scope surfaces or cloud lenses.
- Can be used safely on reusable medical devices, including endoscopes, so there is no need to stock several products.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1623	28 HLD 1 Gallon Bottles	4/case
10-1624	28 HLD 1 Gallon Bottle	1 bottle

GLUTARALDEHYDE STRIPS

PRODUCT NUMBER	FEATURES	QUANTITIES
96-4105	2.1% MEC 60 Strips / Bottle	2/case





SKLAR SOAK™

**INSTRUMENT PRE-SOAK
& DECONTAMINANT
RUST INHIBITOR!**

2 MINUTE KILL TIME

Immersing soiled instruments in SKLAR SOAK effectively pre-cleans and decontaminates, without fear of rusting!

SKLAR SOAK is an effective pre-cleaner that breaks down bio-burden, lifting blood and other organic debris from the surface, while penetrating box locks and crevices. The rust inhibitor protects the delicate chromium oxide layer of the instrument, thus defending against rust and corrosion. Using SKLAR SOAK™ as a pre-cleaner is the first line of defense against rust, prolonging the life of the instrument and therefore protecting your investment!

Soaking for only two minutes thoroughly decontaminates instruments, reducing the risk of exposure to pathogens such as H1N1, HIV-1, HSV, VRE, Rotaviruses, Hep B & C as well as TB, MRSA and A. Baumannii. However, because SKLAR SOAK™ contains a rust inhibitor, instruments may soak for several hours without fear of rusting! Employees are protected from cross contamination, and at the same time can enjoy the freedom to reprocess instruments when it's convenient!

SKLAR SOAK is ideal for use on all stainless steel instruments:

- PH balanced, designed specifically for stainless steel instruments.
- Intermediate level disinfectant with broad spectrum of efficacy.
- Rust inhibitor.
- Non-caustic.
- Ready to use.
- No special handling or disposal.
- Fragrant, no offensive odors.
- Non-glutaraldehyde.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1603	1 Gallon Bottles	4/case
10-1609	1 Gallon Bottle	1 bottle



SKLAR DISINFECTANT

READY-TO-USE SPRAY SURFACE
DISINFECTANT

2 MINUTE KILL TIME

- Ready to use, non-toxic formula.
- Specially formulated for the control of cross contamination.
- Effective against a broad range of pathogens including H1N1, HIV-1, HSV, VRE, Rotaviruses, Hep B & C as well as TB, MRSA and A. Baumannii.
- Safe for all areas of the facility including the NICU, operating rooms, isolation rooms, patient care areas and laboratories.
- Satisfies OSHA's Bloodborne Pathogen Standards.
- Non-irritating formula.
- EPA registered.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1643	24 oz. Spray Bottles	6/case
10-1642	24 oz. Spray Bottle	1 bottle
10-1653	1 Gallon Bottles	4/case
10-1654	1 Gallon Bottle	1 bottle
10-1644	2.5 Gallon Bottles	2/case



SKLAR DISINFECTANT SURFACE WIPES

READY-TO-USE, PRE-MOISTENED
TOWELETTES

Same Effective Formula as Sklar Disinfectant!

2 MINUTE KILL TIME

- Ready to use, non-toxic formula.
- Specially formulated for the control of cross contamination.
- Effective against a broad range of pathogens including H1N1, HIV-1, HSV, VRE, Rotaviruses, Hep B & C as well as TB, MRSA and A. Baumannii.
- Safe for all areas of the facility.
- Satisfies OSHA's Bloodborne Pathogen Standards.
- Destroys odor-causing bacteria.
- 7" x 10" non-woven, tear resistant towelette.
- Clean scent.
- EPA registered.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1617	100 Towelettes / Canister	12/case
10-1616	100 Towelettes / Canister	6/case



SKLAR SPRAY-ZYME™

SURGICAL INSTRUMENT PRE-CLEANER

Multi-Enzyme Based Instrument Presoak & Foam Spray Cleaner

The newest product from Sklar for expert care of your surgical instruments. This multi-enzymatic spray creates a heavy foam which lasts for hours, keeping soiled instruments moist to avoid encrustation while also pre-cleaning prior to processing.

- Prevents Organic Encrustation - enzyme foam stops blood and bodily fluids from drying onto surgical instruments.
- No More Spills - messy spills between the operating room and central cleaning /processing department are eliminated.
- Neutral pH - safe and effective for stainless steel, tungsten carbide, aluminum and brass instruments.
- Bacteriostatic - microbiological control agents inhibit bacterial growth in the ready-to-use foam solution.
- Long Lasting Foam - saves time, foam lasts many hours after application.
- Multi-Enzymatic - contains a synergistic blend of lipase and protease enzymes.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2722	22 oz. Refillable Spray Bottles	12/case
10-2724	1 Gallon Bottles	4/case





SKLAR ENZYMATIC™ DETERGENT

MULTI-ENZYME DETERGENT CONCENTRATE

100% Biodegradable

SKLAR ENZYMATIC Detergent is a low-foaming, fast acting detergent specifically designed for use when instruments and/or cloth materials become soiled with organic debris.

SKLAR ENZYMATIC Detergent may be used as a manual instrument cleaner as well as a pre-soak for blood-contaminated linens, clothing and instruments. It is also used in ultrasonic units and in washer/sterilizer/decontaminators. When used as directed, the enzymatic pH neutral formula is safe for use on rubber, plastic, stainless steel instruments, and surgical equipment as well as cloth.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2777	1 Gallon Bottles	4/case
10-2775	1 Gallon Bottles	2/case
10-2776	1 Gallon Bottle	1 bottle



SKLAR KLEEN™

LOW FOAM LIQUID DETERGENT

The SKLAR KLEEN Low Foam Liquid is concentrated in order to provide a very cost-effective cleaning solution. It is especially formulated to maintain low foam, ensuring effective and trouble-free action in washer/ sterilizers and tunnel washers. Dispersing quickly, the SKLAR KLEEN Low Foam formula is also well suited for automatic injection systems.

- Safe for use in washer / sterilizers.
- Low, Non-Aggressive pH.
- Disperses quickly.
- Non-etching to glassware.
- Concentrated formula.
- Cost-effective solution.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2701	Low Foam 1 Gallon Bottles	4/case
10-2702	Low Foam 1 Gallon Bottle	1 bottle





SKLAR KLEEN™

LIQUID DETERGENT

The SKLAR KLEEN Liquid Detergent was specifically developed for manual and ultrasonic cleaning. It is an effective yet non-corrosive, non aggressive, phosphate free and biodegradable formula that will not etch glassware or attack instrument metals. It is moderately alkaline, affording optimum cleaning power.

- Ideal for manual and ultrasonic cleaning.
- Non-Corrosive.
- Phosphate free and biodegradable.
- Superior wetting power to penetrate locks, serrations, and crevices.
- Appropriate for stainless steel, metal, glassware, plastics, and rubber.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1612	Liquid 8 oz. Bottles	6/case
10-1614	Liquid 8 oz. Bottle	1 bottle
10-1613	Liquid 1 Gallon Bottles	4/case
10-1627	Liquid 1 Gallon Bottle	1 bottle



SKLAR KLEEN™

POWDERED DETERGENT

The SKLAR KLEEN Powdered Detergent is an excellent way to reduce or neutralize the effects of corrosive minerals that are inherent in virtually all community water supplies. SKLAR KLEEN Powdered Detergent is ideal for manual and ultrasonic cleaning. One economical 3½ lb. container can make in excess of 150 gallons of cleaning solution. It's non-corrosive, non aggressive, phosphate free formula will not etch glassware or attack your instrument metals.

- Economical and Convenient.
- For all manual and ultrasonic cleaning.
- Appropriate for stainless steel, metal glassware, plastics and rubber.
- Non-Corrosive.
- Easy to Handle.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1608	Powdered 5 lb. Containers	4/case
10-1607	Powdered 5 lb. Container	1 container
10-1630	Powdered 3.5 lb. Pitchers	6/case
10-2765	Powdered 3.5 lb. Pitcher	1 pitcher





SKLAR INSTRU-GUARD™ ONE STEP

**A "ONE STEP" INSTRUMENT
CLEANER AND LUBRICANT**

EXTENDS INSTRUMENT LIFE!

SKLAR INSTRU-GUARD ONE STEP contains a penetrating lubricant that protects and helps extend the life of surgical instruments. The detergents and emulsifiers quickly break up and remove organic soils from box locks and other hard to reach places.

SKLAR INSTRU-GUARD ONE STEP contains four rust inhibitors to control corrosion and eliminate scaling and spotting. It does not contain phosphates and can reduce the time and cost involved in processing your instruments by cleaning and lubricating in one easy step.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2800	1 Gallon Bottles	4/case
10-2801	1 Gallon Bottle	1 bottle



SKLAR INSTRU-GUARD™ LUBE

SURGICAL INSTRUMENT LUBRICANT / RUST INHIBITOR

- Increase the Life of Instruments.
- Protects all Surgical Instruments.

SKLAR INSTRU-GUARD LUBE is an easy-to-use, water soluble, lubricant and rust inhibitor. When used daily as a rinse after every cleaning it keeps moving parts operating smoothly and remains effective after autoclaving, but without leaving an oily residue.

Regular lubrication is essential no matter which cleaning method is employed (including ultrasonic). SKLAR INSTRU-GUARD LUBE is ideal for maintaining instrument quality and performance.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1635	1 Gallon Bottles Concentrated	4/case
10-1634	1 Gallon Bottle Concentrated	1 bottle
10-1636	8 oz. Spray Bottles	12/case
10-1637	8 oz. Spray Bottle	1 bottle





SKLAR POLISH™

INSTRUMENT STAIN REMOVER

SKLAR POLISH is an effective stain remover for use on all your metal surgical instruments, including stainless steel, chrome and brass. It will brighten instrument trays, basins and racks and may also be used on porcelain and tile surfaces.

- Non Abrasive.
- Removes corrosive stains.
- Extends instruments life.
- "Like New" results.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1626	8 oz. Jars	12/case
10-1629	8 oz. Jars	6/case
10-1927	8 oz. Jar	1 jar



INSTRUMENT CLEANING BRUSHES

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1650	Stainless Steel Bristles	3/pk
10-1652	Stainless Steel Bristles	12/pk
10-1658	Stainless Steel Bristles	24/pk
10-1665	Stainless Steel Bristles	48/pk
10-1651	Nylon Bristles	3/pk
10-1657	Nylon Bristles	12/pk
10-1659	Nylon Bristles	24/pk
10-1656	Nylon Bristles	48/pk



SKLAR-SHEEN™

AUTOCLAVE AND INSTRUMENT CLEANER AND CONDITIONER

With an efficacious phosphoric acid cleaning agent, Sklar Sheen is a robust spot, stain and rust remover. Ideal for dissolving hard water deposits, use in autoclaves to quickly cut through steam sterilizer scale and deposits.

Versatile and effective, SKLAR SHEEN couples powerful stain and rust removers with corrosion inhibitors, making it tough on stains but gentle on instruments. Recondition instruments to restore them to nearly brand new working condition.

Use to clean and recondition any stainless steel surface to its original luster, including carts, trays and basins.

- Conditions for lasting protection.
- Dissolves hard water deposits.
- Frees up frozen box locks.
- Restores and extends surgical instrument life.

PRODUCT NUMBER	FEATURES	QUANTITIES
10-2738	22 oz. Spray Bottles	12/case
10-2736	22 oz. Spray Bottle	1 bottle



SKLAR INSTRUMENT CARE KIT

A TRIO OF SKLAR'S CARE AND CLEANING PRODUCTS
TO KEEP INSTRUMENTS IN TIPTOP SHAPE

PRODUCT NUMBER	FEATURES	QUANTITIES
10-1625		
	Sklar Lube™ 8 oz. Spray Bottle	1 each
	Sklar Kleen™ Liquid 8 oz. Bottle	1 each
	Sklar Polish™ 8 oz. Jar	1 each
	Stainless Steel Bristle Cleaning Brush	1 each
	Nylon Bristle Cleaning Brush	1 each



TROUBLE SHOOTING

NO STEEL IS TRULY “STAINLESS”

It's a fact that so-called stainless steel is subject to both water-spotting and staining. However what many people identify as a “rust” problem usually turns out to be something else.

Most often the problem is actually a stain caused by a surface deposit. These deposits appear in a variety of colors depending upon the type of deposit involved. Once the cause is found and eliminated, most of the time the problem quickly diminishes or disappears completely.

As you probably know, stainless steel is ideally suited for the surgical suite because it is rust resistant, it can

be honed to an extremely sharp edge or fine point, and it can be hardened to maintain the delicate, yet precise, requirements of the surgeon. But it is this hardenable quality – the amount of carbon used during the manufacturing process – coupled with harsh chemicals and improper care that renders stainless steel susceptible to corrosion.

To protect your instrument investment, Sklar cleaning agents have been specifically formulated to provide cost effective instrument care solutions that help maintain the life and performance of your valuable stainless steel instruments.



889 South Matlack Street

West Chester, PA 19382

610.430.3200

Fax: 610.696.9007

service@sklarcorp.com

www.sklarcorp.com

©Copyright 2016. Sklar Corporation, West Chester, PA. All rights reserved. No portion of this catalog may be reproduced without written permission of Sklar.